

1


Assemble Sizing Chart First

2010

Important

Before printing, turn off "fit to page" in your print options dialogue box.*

A



B

Fold under on *Solid Fold Line* on page 2, align *Solid Lines* and arrows, and tape pages together.


2

Length Measurement

A

- Fold sizing chart on dotted line
- Make a 90-degree angle
- Place on hard surface
- Put fold firmly against wall
- Tape to prevent mis-sizing

use arrows to align pages



Children to adult sizes

use arrows to align pages

Solid Line

use
arrows
to align
pages

B

- Place heel of AFO on sizing chart
- Be certain to include posterior stop
- If foot is longer than brace, include foot

C

- Trace AFOs for accuracy
- Match longest point to number on right
- If you are between or slightly above line, round up


D

IMPORTANT
Necessary Adjustments for Circumference

For all AFOs, TCFlex, DAFOs, KAFOs, SMOs with:

- Solid ankle, instep strapping
Increase 1 full size
- Articulating hinges, posterior stop and instep strapping
Increase 1.5
- Turbo boot, posterior stop, hinges, footplates, instep, forefoot and toe strapping;
Increase 2 sizes

Fold under on *Solid Fold Line*, align *Solid Lines* and arrows, and tape pages together.



Fold up on *Dotted Fold Line*.

KEEPING PACE

Children's Orthopedic Footwear



OrtoPed
www.ortoped.ca
373 McCaffrey St.
Montréal, QC
Canada, H4T 1Z7
Tel: 514.342.6161
or 1.800.363.8726
Fax: 514.342.7565
or 1.800.663.8817
info@ortoped.ca
www.ortoped.ca

use
arrows
to align
pages