

» General Purpose Gel Liner

Liner

Supplier:

ALPS South Corporation

Features & Benefits:

- High performance fabric improves durability
- Reduced bunching in the popliteal region
- Re-engineered distal end encapsulation can retrofit most liners
- Easier donning with "slick fabric"

Specifications:

- Types: *Locking & Cushion liners*
- Style: *Uniform*
- Material: *ALPS EasyGel*
- Thicknesses: *3mm & 6mm*
- Sizes: *Eight (See chart)*
- Measurements taken 6cm from distal end

Indications:

- Residual limbs needing increased cushioning and shear relief

Contraindications:


- Poor hand dexterity and strength for donning

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- 12 months for single liner or 24 months for two identical liners, against manufacturers defects



General Purpose Liner » ALPS South Corporation

Item No.	Description	Style	Thickness	Size	Distal Circumference*	Length
ALPS-GPDT-16-3	Locking gel liner	Uniform	3mm	16	16–19 cm	13–14"
ALPS-GPDT-20-3	Locking gel liner	Uniform	3mm	20	20–23 cm	17–18"
ALPS-GPDT-24-3	Locking gel liner	Uniform	3mm	24	24–25 cm	17–18"
ALPS-GPDT-26-3	Locking gel liner	Uniform	3mm	26	26–27 cm	17–18"
ALPS-GPDT-28-3	Locking gel liner	Uniform	3mm	28	28–31 cm	17–18"
ALPS-GPDT-32-3	Locking gel liner	Uniform	3mm	32	32–37 cm	17–18"
ALPS-GPDT-38-3	Locking gel liner	Uniform	3mm	38	38–43 cm	17–18"
ALPS-GPDT-44-3	Locking gel liner	Uniform	3mm	44	44–53 cm	17–18"
ALPS-GPDT-16-6	Locking gel liner	Uniform	6mm	16	16–19 cm	13–14"
ALPS-GPDT-20-6	Locking gel liner	Uniform	6mm	20	20–23 cm	17–18"
ALPS-GPDT-24-6	Locking gel liner	Uniform	6mm	24	24–25 cm	17–18"
ALPS-GPDT-26-6	Locking gel liner	Uniform	6mm	26	26–27 cm	17–18"
ALPS-GPDT-28-6	Locking gel liner	Uniform	6mm	28	28–31 cm	17–18"
ALPS-GPDT-32-6	Locking gel liner	Uniform	6mm	32	32–37 cm	17–18"
ALPS-GPDT-38-6	Locking gel liner	Uniform	6mm	38	38–43 cm	17–18"
ALPS-GPDT-44-6	Locking gel liner	Uniform	6mm	44	44–53 cm	17–18"
ALPS-GPFR-16-3	Cushion gel liner	Uniform	3mm	16	16–19 cm	13–14"
ALPS-GPFR-20-3	Cushion gel liner	Uniform	3mm	20	20–23 cm	17–18"
ALPS-GPFR-24-3	Cushion gel liner	Uniform	3mm	24	24–25 cm	17–18"
ALPS-GPFR-26-3	Cushion gel liner	Uniform	3mm	26	26–27 cm	17–18"
ALPS-GPFR-28-3	Cushion gel liner	Uniform	3mm	28	28–31 cm	17–18"
ALPS-GPFR-32-3	Cushion gel liner	Uniform	3mm	32	32–37 cm	17–18"
ALPS-GPFR-38-3	Cushion gel liner	Uniform	3mm	38	38–43 cm	17–18"
ALPS-GPFR-44-3	Cushion gel liner	Uniform	3mm	44	44–53 cm	17–18"
ALPS-GPFR-16-6	Cushion gel liner	Uniform	6mm	16	16–19 cm	13–14"
ALPS-GPFR-20-6	Cushion gel liner	Uniform	6mm	20	20–23 cm	17–18"
ALPS-GPFR-24-6	Cushion gel liner	Uniform	6mm	24	24–25 cm	17–18"
ALPS-GPFR-26-6	Cushion gel liner	Uniform	6mm	26	26–27 cm	17–18"
ALPS-GPFR-28-6	Cushion gel liner	Uniform	6mm	28	28–31 cm	17–18"
ALPS-GPFR-32-6	Cushion gel liner	Uniform	6mm	32	32–37 cm	17–18"
ALPS-GPFR-38-6	Cushion gel liner	Uniform	6mm	38	38–43 cm	17–18"
ALPS-GPFR-44-6	Cushion gel liner	Uniform	6mm	44	44–53 cm	17–18"

* Measure distal end circumference approximately 6cm from distal stump end