

## » Thermoliner® Custom Liner

### Liner

#### Supplier:

ALPS South Corporation

#### Features & Benefits:

- Can be custom molded at prosthetist facility. No cast to send.
- Prosthetist can customize and shape the liner to the configuration of the user
- Can be pulled on like a sock, making donning easier
- Gel durometer and contour provide additional control and support
- Provides solution for difficult-to-fit patients
- Adaptable for most types of users

#### Specifications:

- Types: *Cushion and Locking Liners*
- Style: *Uniform*
- Material: *Thermoformable gel*
- Thicknesses: *3mm and 6mm*
- Sizes: *16, 20, 24, 26, 28, 32, 38 and 44*
- Locking liners come with distal attachment
- Measurements taken 6cm from the distal end

#### Indications:


- Transtibial and transfemoral amputations
- Symes/Knee disarticulation amputation
- Athletes
- Difficult shapes and sizes

#### Call Points:

- Prosthetist
- Prosthetic technician

#### Warranty:

- 13 months, against manufacturers defects



**Thermoliner® Custom Liner » ALPS South Corporation**

Item No.	Description	Style	Thickness	Size	Distal Circumference*	Length
<b>Cushion Liner</b>						
ALPS-TFFR-16-3	Thermoformable cushion gel liner	Uniform	3mm	16	16–19 cm	13–14"
ALPS-TFFR-20-3	Thermoformable cushion gel liner	Uniform	3mm	20	20–23 cm	17–18"
ALPS-TFFR-24-3	Thermoformable cushion gel liner	Uniform	3mm	24	24–25 cm	17–18"
ALPS-TFFR-26-3	Thermoformable cushion gel liner	Uniform	3mm	26	26–27 cm	17–18"
ALPS-TFFR-28-3	Thermoformable cushion gel liner	Uniform	3mm	28	28–31 cm	17–18"
ALPS-TFFR-32-3	Thermoformable cushion gel liner	Uniform	3mm	32	32–37 cm	17–18"
ALPS-TFFR-38-3	Thermoformable cushion gel liner	Uniform	3mm	38	38–43 cm	17–18"
ALPS-TFFR-44-3	Thermoformable cushion gel liner	Uniform	3mm	44	44–53 cm	17–18"
ALPS-TFFR-16-6	Thermoformable cushion gel liner	Uniform	6mm	16	16–19 cm	13–14"
ALPS-TFFR-20-6	Thermoformable cushion gel liner	Uniform	6mm	20	20–23 cm	17–18"
ALPS-TFFR-24-6	Thermoformable cushion gel liner	Uniform	6mm	24	24–25 cm	17–18"
ALPS-TFFR-26-6	Thermoformable cushion gel liner	Uniform	6mm	26	26–27 cm	17–18"
ALPS-TFFR-28-6	Thermoformable cushion gel liner	Uniform	6mm	28	28–31 cm	17–18"
ALPS-TFFR-32-6	Thermoformable cushion gel liner	Uniform	6mm	32	32–37 cm	17–18"
ALPS-TFFR-38-6	Thermoformable cushion gel liner	Uniform	6mm	38	38–43 cm	17–18"
ALPS-TFFR-44-6	Thermoformable cushion gel liner	Uniform	6mm	44	44–53 cm	17–18"
<b>Locking Liner</b>						
ALPS-TFDT-16-3	Thermoformable locking gel liner	Uniform	3mm	16	16–19 cm	13–14"
ALPS-TFDT-20-3	Thermoformable locking gel liner	Uniform	3mm	20	20–23 cm	17–18"
ALPS-TFDT-24-3	Thermoformable locking gel liner	Uniform	3mm	24	24–25 cm	17–18"
ALPS-TFDT-26-3	Thermoformable locking gel liner	Uniform	3mm	26	26–27 cm	17–18"
ALPS-TFDT-28-3	Thermoformable locking gel liner	Uniform	3mm	28	28–31 cm	17–18"
ALPS-TFDT-32-3	Thermoformable locking gel liner	Uniform	3mm	32	32–37 cm	17–18"
ALPS-TFDT-38-3	Thermoformable locking gel liner	Uniform	3mm	38	38–43 cm	17–18"
ALPS-TFDT-44-3	Thermoformable locking gel liner	Uniform	3mm	44	44–53 cm	17–18"
ALPS-TFDT-16-6	Thermoformable locking gel liner	Uniform	6mm	16	16–19 cm	13–14"
ALPS-TFDT-20-6	Thermoformable locking gel liner	Uniform	6mm	20	20–23 cm	17–18"
ALPS-TFDT-24-6	Thermoformable locking gel liner	Uniform	6mm	24	24–25 cm	17–18"
ALPS-TFDT-26-6	Thermoformable locking gel liner	Uniform	6mm	26	26–27 cm	17–18"
ALPS-TFDT-28-6	Thermoformable locking gel liner	Uniform	6mm	28	28–31 cm	17–18"
ALPS-TFDT-32-6	Thermoformable locking gel liner	Uniform	6mm	32	32–37 cm	17–18"
ALPS-TFDT-38-6	Thermoformable locking gel liner	Uniform	6mm	38	38–43 cm	17–18"
ALPS-TFDT-44-6	Thermoformable locking gel liner	Uniform	6mm	44	44–53 cm	17–18"

\* Measure distal end circumference approximately 6cm from distal stump end