

» Combo

Dynamic KAFO

Supplier:

Allard Inc.

Features & Benefits:

- Low profile and lightweight
- Uprights easily shaped and adjusted to optimum height for patient
- Modular components allow you to accommodate variances in thigh and calf circumferences
- Easy to don and doff when patient is resting
- Can be fabricated onto ToeOFF® or BlueRocker™ within an hour

Specifications:

- Weight limit: 150kg (330lbs)
- Height limit: 2m (6'6")
- Tibia connector sizes: *Small, Medium and Large/X-Large (see chart)*
- Knee control component sizes: *Small, Medium and Large (see chart)*
- Thigh cuff sizes: *Small, Medium, Large and X-Large (see chart)*
- For use with ToeOFF® or BlueRocker™ only
- Fitting requires a trained professional
- Ordering of this item requires a selection of the following components to be made: Joint and Upright, Tibia Connector, Knee Control Component, Thigh Cuff and SoftKIT™

Indications:

- Management of genu recurvatum or chronic knee instability, accompanied with footdrop

Contraindications:

- Total knee ligament loss
- Complete quadriceps insufficiency
- Severe loss of cognitive skills

Call Points:

- Orthotist
- Orthotic technician

Warranty:

- One year from date of initial fitting on ToeOFF® or BlueRocker™, against manufacturer defects



Combo » Allard Inc.

Please note: When ordering for this item, please make a selection of the following components listed below:

Item No.	Description	Size
Joints with Uprights		
AL-287800000	Dynamic KAFO joints and uprights	One size

Item No.	Description	Size
Tibia Connector		
AL-287840011	Dynamic KAFO tibia connector	Small
AL-287840012	Dynamic KAFO tibia connector	Medium
AL-287840013	Dynamic KAFO tibia connector	Large/X-Large

Please note: Generally, select the same size as the ToeOFF or BlueRocker. However, the next smaller size will "snug" the wings of the device closer to the leg for the extra small calf, and the next larger size may better accommodate a larger calf.

Item No.	Description	Size	Medial /Lateral Dia. ¹
Knee Control Component with popliteal interface			
AL-287820011	Dynamic KAFO knee control component	Small	70–105 mm (2¾"–4¼")
AL-287820012	Dynamic KAFO knee control component	Medium	95–130 mm (3¾"–5¼")
AL-287820013	Dynamic KAFO knee control component	Large	125–160 mm (5"–6¼")

¹For proper fitting, take medial/lateral diameter measurement at the knee

Item No.	Item No.	Description	Size	Circumference ²	Medial /Lateral Dia. ²
Left	Right				
Thigh Cuff					
AL-28711011	AL-28712011	Dynamic KAFO thigh cuff	Small	30–40 cm (11¾"–15¾")	11–15 cm (4¼"–6")
AL-28711012	AL-28712012	Dynamic KAFO thigh cuff	Medium	38–50 cm (15–19¾")	13–17 cm (5¼"–6¾")
AL-28711013	AL-28712013	Dynamic KAFO thigh cuff	Large	48–60 cm (19–23¾")	15–19 cm (6"–7½")
AL-28711014	AL-28712014	Dynamic KAFO thigh cuff	X-Large	58–70 cm (22¾"–27½")	17–21 cm (6¾"–8¼")

²For proper fitting, take circumference and medial/lateral diameter measurements 10cm (4") above tibial plateau

Combo » Padding

Item No.	Description	Size	Color
AL-287830011	SoftKIT™ for Combo Dynamic KAFO thigh cuff — Set of two	Small	Black
AL-287830012	SoftKIT™ for Combo Dynamic KAFO thigh cuff — Set of two	Medium	Black
AL-287830013	SoftKIT™ for Combo Dynamic KAFO thigh cuff — Set of two	Large	Black
AL-287830014	SoftKIT™ for Combo Dynamic KAFO thigh cuff — Set of two	X-Large	Black