

» HeelWedge™

Medical Shoe

Supplier:

Darco International Inc.

Features & Benefits:

- Square toe design provides additional room to protect toes
- Wide base for extra stability
- Removable forefoot closure provides security while eliminating buckle pressure
- Ankle strap seats foot firmly in shoe to reduce heel slippage and eliminate friction
- Zoned outsole has an aggressive tread under the midfoot where traction is needed most
- Removable insole features twice the padding of standard insoles and can be modified or replaced with the PegAssist™ off-loading insole

Specifications:

- Sizes: *X-Small, Small, Medium, Large, and X-Large*
- Color: *Black*
- It is suggested that the period of time for patient use is 8 to 12 weeks

Indications:

- Rearfoot trauma
- Wounds or ulcerations present on the heel area
- Post surgical healing for either soft tissue or bony structure of the heel

Call Points:

- Podiatrist
- Pedorthist
- Orthotist

Warranty:

- Against manufacturer defects



Item No.	Description	Size	Size (Mens)*	Size (Womens)*	Color
DAR-HQ0B	Rearfoot off-loading medical shoe	X-Small	—	4–7	Black
DAR-HQ1B	Rearfoot off-loading medical shoe	Small	6–8	7½–10	Black
DAR-HQ2B	Rearfoot off-loading medical shoe	Medium	8½–10	10½–13	Black
DAR-HQ3B	Rearfoot off-loading medical shoe	Large	10½–12	13+	Black
DAR-HQ4B	Rearfoot off-loading medical shoe	X-Large	12½–14	—	Black

* Approximate sizes only. For proper fitting, consult with patient or measure their shoe size prior to ordering.