

**HeelWedge™
Medical Shoe**

Vendor:
Darco International Inc.

Applications:

- Off-loading
- Woundcare

Patient Benefits:

- Secure
- Safe
- Hygienic

Features:

- Square toe design provides more roomier toe box
- Off-loading for heel area
- Rocker bottom
- Unique forefoot closure without the pressure points of buckles and straps
- Cloth toe protector keeps toes warm and dressing clean and dry
- Ankle strap seats the foot firmly in the sole and eliminates heel friction

Indications:

- Plantar fasciitis
- Ulcerations
- Infection
- Trauma
- Promote healing after surgery or when wounds are present at the heel
- Suggested use: For a period of 8 to 12 weeks

Call Points:

- Podiatrist
- Pedorthist
- Orthotist

Warranty:

- Against manufacturer defects



Darco International Inc., HeelWedge™ Shoe

Item No	Description	Size	Color	Women's Size	Men's Size
DAR-HQ0B	HeelWedge™ Shoe	X-Small	Black	4-7	—
DAR-HQ1B	HeelWedge™ Shoe	Small	Black	7.5-10	6-8
DAR-HQ2B	HeelWedge™ Shoe	Medium	Black	10-13	8.5-11
DAR-HQ3B	HeelWedge™ Shoe	Large	Black	—	10.5-12
DAR-HQ4B	HeelWedge™ Shoe	X-Large	Black	—	12.5-14