

## » OrthoWedge™

### Medical Shoe

#### Supplier:

Darco International Inc.

#### Features & Benefits:

- Square toe design provides additional room to protect toes
- Wide base for extra stability
- Ankle strap seats foot firmly in shoe to reduce heel slippage and eliminate friction
- Removable forefoot closure provides security while eliminating buckle pressure
- Removable insole features twice the padding of standard insoles and can be modified or replaced with the PegAssist™ off-loading insole

#### Specifications:

- Sizes: *X-Small, Small, Medium, Large, and X-Large*
- Color: *Black*
- It is suggested that the period of time for patient use is 8 to 12 weeks

#### Indications:

- Post surgical healing
- Forefoot trauma
- Wounds or ulcerations present under metatarsal heads and toes
- Situations that requires off-loading pressure of the forefoot

#### Contraindications:

- Constant excessive weight
- Walking on uneven terrain

#### Call Points:

- Podiatrist
- Pedorthist
- Orthotist

#### Warranty:

- Against manufacturer defects



Item No.	Description	Size	Size (Mens)*	Size (Womens)*	Color
DAR-0Q0B	Forefoot off-loading medical shoe	X-Small	—	4–7	Black
DAR-0Q1B	Forefoot off-loading medical shoe	Small	6–8	7½–10	Black
DAR-0Q2B	Forefoot off-loading medical shoe	Medium	8½–10	10½–13	Black
DAR-0Q3B	Forefoot off-loading medical shoe	Large	10½–12	13+	Black
DAR-0Q4B	Forefoot off-loading medical shoe	X-Large	12½–14	—	Black

\* Approximate sizes only. For proper fitting, consult with patient or measure their shoe size prior to ordering.