

OrthoWedge™
Medical Shoe

Vendor:

Darco International Inc.

Applications:

- Reduce pressure from forefoot, woundcare, post-op

Patient Benefits:

- Stability

Features:

- Square toe design provides additional room to protect toes
- Wide base for extra stability
- Ankle strap to eliminate heel slippage
- Removable forefoot closure to eliminate buckle pressure
- 10mm insole, fully customizable with over twice the padding
- 15° dorsi-angle shifts weight back to the midfoot and heel

Indications:

- Any condition from the metatarsal heads distally in which it is desirable to reduce body weight such as ulcers, infections, trauma and surgery
- Suggested use: For a period of 8 to 12 weeks

Contraindications:

- Constant excessive weight
- Walking on uneven terrain

Call Points:

- Podiatrist
- Pedorthist
- Orthotist

Warranty:

- Against manufacturer defects



Darco International Inc., OrthoWedge™ Shoe

Item No.	Description	Size	Color	Women's Sizes	Men's Sizes
DAR-OQ0B	OrthoWedge™ Medical Shoe	X-Small	Black	4–7	—
DAR-OQ1B	OrthoWedge™ Medical Shoe	Small	Black	7.5–10	6–8
DAR-OQ2B	OrthoWedge™ Medical Shoe	Medium	Black	—	8.5–10
DAR-OQ3B	OrthoWedge™ Medical Shoe	Large	Black	—	10.5–12
DAR-OQ4B	OrthoWedge™ Medical Shoe	X-Large	Black	—	12.5–14