

Peg-Assist™ System

Insole

Supplier:

Darco International Inc.

Features & Benefits:

- Pegs can be removed to off-load pressure from any area of the plantar surface of the foot
- Provides the perfect combination of shock absorption and moldability
- After peg removal, the Poron® cover remains intact to eliminate the incidence of ring edema and assist in supporting the area
- Removable toe cover helps to ensure that the forefoot remains clean, dry and free of foreign objects
- Perfect choice for the treatment of diabetic ulcerations
- Provide short-term pain management for patients suffering due to bony abnormalities (e.g. prominent metatarsal heads) or other painful conditions that may later be
- Eliminate “lever effect” on osteotomies of the MT1 by off-loading the first MPJ. The insole can also be cut flush with the metatarsal heads to “float” the toes when K-Wires are present
- A single-use marking tube is included to mark the area of desired pressure reduction on patients foot and to ensure proper peg removal

Specifications:

- Thickness: 18 mm
- Combination of 5 mm Plastazote®, 3 mm Poron® and 10 mm multi-foam
- Available in 5 different sizes (see chart)

Indications:

- Off-loading
- Wound care
- Post-op care
- Pain management
- Foot trauma

Call Points:

- Orthotist
- Physician
- Orthopedic surgeon
- Cast technician
- Podiatrist

Warranty:

- Against manufacturer defects



Peg-Assist™ System » Darco International Inc.

Item No.	Description	Gender	Size	Shoe Size	Color
<i>MEN</i>					
DAR-PTQM2	Foot pressure off-loading system	Men	Medium	8.5–10	Black/Blue/Pink
DAR-PTQM3	Foot pressure off-loading system	Men	Large	10.5–12	Black/Blue/Pink
DAR-PTQM4	Foot pressure off-loading system	Men	X-Large	12.5–14	Black/Blue/Pink
<i>WOMEN</i>					
DAR-PTQW1	Foot pressure off-loading system	Women	Small	4–6	Black/Blue/Pink
DAR-PTQW2	Foot pressure off-loading system	Women	Medium	6.5–8	Black/Blue/Pink
DAR-PTQW3	Foot pressure off-loading system	Women	Large	8.5–10	Black/Blue/Pink