

» Wound Care Shoe System™

Medical Shoe

Supplier:

Darco International Inc.

Features & Benefits:

- Off-loads areas of the foot while patient remains ambulatory
- Four multi-density insoles allow for customization in precisely targeted off-loading
- Leather upper can be modified or cut away to provide even greater levels of pressure relief
- Deep rocker sole for easy ambulation
- Easy hook and loop closures

Specifications:

- Sizes: *X-Small, Small, Small-Medium, Medium, Large, X-Large and XX-Large*
- Materials: *Shoe= Leather with Plastazote® lining*
Insoles= Various densities of material with Plastazote® top
- Color: *Black*
- The suggested period of time for patient use is 8 to 12 weeks

Indications:


- Diabetic foot ulcers and infections

Call Points:

- Podiatrist
- Pedorthist
- Orthotist

Warranty:

- Against manufacturer defects



Insoles for the Woundcare Shoe System

Item No.	Description	Size	Size (Mens)*	Size (Womens)*	Color
DAR-WCS0B	Woundcare shoe	X-Small	—	4–5½	Black
DAR-WCS1B	Woundcare shoe	Small	6–7½	6–8½	Black
DAR-WCS7B	Woundcare shoe	Small-Medium	8–8½	9–9½	Black
DAR-WCS2B	Woundcare shoe	Medium	9–9½	10–10½	Black
DAR-WCS3B	Woundcare shoe	Large	10–11½	11–12½	Black
DAR-WCS4B	Woundcare shoe	X-Large	12–13½	13–14½	Black
DAR-WCS5B	Woundcare shoe	XX-Large	14–16	15–17	Black

* Approximate sizes only. For proper fitting, consult with patient or measure their shoe size prior to ordering.