

**Wound Care System**  
**Medical Shoe**

**Vendor:**

Darco International Inc.

**Practitioner Benefits:**

- Easily off-load areas of the foot while the patient remains ambulatory

**Features:**

- Easy hook and loop closures
- Deep rocker sole to provide easy ambulation
- Sole contains 4 layers of differing density insoles that may be altered for off-loading
- Leather upper

**Specifications:**

- Leather upper is lined with plastazote® material

**Indications:**

- Permits off-loading ulcers, infected area, bony prominence on both the plantar surfaces of the foot
- Diabetic foot ulcers and infections
- Suggested use: For a period of 8 to 12 weeks

**Call Points:**

- Podiatrist
- Pedorthist
- Orthotist

**Warranty:**

- Against manufacturer defects



**Darco International Inc., Woundcare Shoe System**

Item No.	Description	Size	Color	Woman's Size	Men's Size
DAR-WCS0B	Woundcare Shoe	X-Small	Black	4-5.5	—
DAR-WCS1B	Woundcare Shoe	Small	Black	6-8.5	6-7.5
DAR-WCS7B	Woundcare Shoe	Sm-Medium	Black	9-9.5	8-8.5
DAR-WCS2B	Woundcare Shoe	Medium	Black	10-10.5	9-9.5
DAR-WCS3B	Woundcare Shoe	Large	Black	11-12.5	10-11.5
DAR-WCS4B	Woundcare Shoe	X-Large	Black	13-14.5	12-13.5
DAR-WCS5B	Woundcare Shoe	XX-Large	Black	15-17	14-16