

» EZ-Crawl

Talipes Equinovarus Orthosis

Supplier:

Fillauer LLC

Features & Benefits:

- Treatment of congenital clubfoot, using Ponsetti Method
- Designed to be fixated to custom AFO's, above malleoli
- Quick disconnects on bar allows AFO's to be removed or attached quickly for diaper changes, donning or doffing
- Permits standing and full weight bearing
- Adjustable bar width to accommodate growth
- Allows footwear with orthosis
- Lightweight and versatile design

Specifications:

- Material: *Bar= Aluminum*
Anchor points= Stainless steel
- Minimum length: 254mm (10")
- Maximum length: 406mm (16")
- Width: 25mm (1")
- Includes mounting kit and Jungle Animal II transfer paper (DF-999102)
- Involved foot (or feet) should be positioned with 70° of external rotation. The unaffected foot should be maintained at 40° of external rotation.

Indications:

- Congenital Talipes Equinovarus (clubfoot)

Call Points:

- Orthotist
- Orthotic Technician

Warranty:

- Against manufacturers defects



Jungle Animal II transfer paper (DF-999102) included

Item No.	Description	Min. Length	Max. Length	Width
DF-012400	Talipes equinovarus orthosis	254mm (10")	406mm (16")	25mm (1")