

» **Ibex / Ibex XD Pylon Foot System™**

Prosthetic Foot

Supplier:

Fillauer LLC

Features & Benefits:

- Unique Micro-Slice technology enhances inversion and eversion of the forefoot
- Full length angled upper spring provides fluid gait, increased energy storage and dynamic response
- Interchangeable heel elastomers for fine tuning heel stiffness
- Modular design with integrated pyramid
- Durable and lightweight
- Ibex XD allows for greater weight bearing at 200kg
- Innovative foot shell — flexible with realistic features

Specifications:

- Models: *Ibex*
Ibex XD (high load bearing)
- Sizes: *Ibex* = 22–30 cm
Ibex XD = 25–30 cm
- Weight limit: *Ibex* = 150kg (330lbs)
Ibex XD = 200kg (440lbs)
- Weight: *Ibex* = 459g (16.2oz)*
(Based on size 27cm foot, without foot shell)
Ibex XD = 602g (21.2oz)
(Based on size 30cm foot, without foot shell)
- Spring Stiffnesses: 2–9 (see chart)
- Heel height: 10mm ($\frac{3}{8}$ ")
- Clearance: *Ibex* = 170mm ($6\frac{3}{4}$ ")
Ibex XD = 194mm ($7\frac{5}{8}$ ")
- Foot shell: *Caucasian (custom colors are available for order)*
- Activity level: *High (K3)*
- Complies with ISO-22675 test standard
- Includes heel wedges, Spectra sock and instruction manual

Indications:

- Transfemoral and transtibial amputees
- Active K3 amputees

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- 36 months on foot, against manufacturers defects
- Six months on foot shell, against manufacturers defects



Ibex Pylon Foot System™



Ibex XD Foot System™

Ibex / Ibex XD Pylon Foot System™ » Fillauer LLC

Ordering example: DF-120-10-25-E8XD

Item No.	Description	Size	Clearance	Heel Height	Weight ¹
DF-120-10-22- KK	High impact prosthetic foot	22cm	170mm (6¾")	10mm (¾")	459g (16.2oz) ¹
DF-120-10-23- KK	High impact prosthetic foot	23cm	170mm (6¾")	10mm (¾")	459g (16.2oz) ¹
DF-120-10-24- KK	High impact prosthetic foot	24cm	170mm (6¾")	10mm (¾")	459g (16.2oz) ¹
DF-120-10-25- KK	High impact prosthetic foot	25cm	170mm (6¾")	10mm (¾")	459g (16.2oz) ¹
DF-120-10-26- KK	High impact prosthetic foot	26cm	170mm (6¾")	10mm (¾")	459g (16.2oz) ¹
DF-120-10-27- KK	High impact prosthetic foot	27cm	170mm (6¾")	10mm (¾")	459g (16.2oz) ¹
DF-120-10-28- KK	High impact prosthetic foot	28cm	170mm (6¾")	10mm (¾")	459g (16.2oz) ¹
DF-120-10-29- KK	High impact prosthetic foot	29cm	170mm (6¾")	10mm (¾")	459g (16.2oz) ¹
DF-120-10-30- KK	High impact prosthetic foot	30cm	170mm (6¾")	10mm (¾")	459g (16.2oz) ¹
DF-120-10-25- KKXD	High load bearing, high impact prosthetic foot	25cm	194mm (7⅝")	10mm (¾")	602g (21.2oz) ²
DF-120-10-26- KKXD	High load bearing, high impact prosthetic foot	26cm	194mm (7⅝")	10mm (¾")	602g (21.2oz) ²
DF-120-10-27- KKXD	High load bearing, high impact prosthetic foot	27cm	194mm (7⅝")	10mm (¾")	602g (21.2oz) ²
DF-120-10-28- KKXD	High load bearing, high impact prosthetic foot	28cm	194mm (7⅝")	10mm (¾")	602g (21.2oz) ²
DF-120-10-29- KKXD	High load bearing, high impact prosthetic foot	29cm	194mm (7⅝")	10mm (¾")	602g (21.2oz) ²
DF-120-10-30- KKXD	High load bearing, high impact prosthetic foot	30cm	194mm (7⅝")	10mm (¾")	602g (21.2oz) ²

KK = Spring. See the below chart for spring selection.

¹Based on a 26cm foot, without foot shell

²Based on a 30cm foot, without foot shell

Ibex / Ibex XD » Spring Stiffness (KK)

	45–54 kg (100–119 lbs)	54–63 kg (120–139 lbs)	64–72 kg (140–159 lbs)	73–81 kg (160–179 lbs)	82–90 kg (180–199 lbs)	91–102 kg (200–224 lbs)	102–113 kg (225–249 lbs)	113–125 kg (250–275 lbs)	125–136 kg (276–300 lbs)	137–150 kg (301–330 lbs)	150–200 kg (331–440 lbs)
22cm	A2	A2	B3	B3	B4	B4	—	—	—	—	—
23cm	A2	A2	B3	B3	B4	B4	—	—	—	—	—
24cm	A2	A2	B3	B3	B4	B4	—	—	—	—	—
25cm	B4	B4	B5	B5	C6	C6	D7	D7	E8	E8	E8 ¹
26cm	B4	B4	B5	B5	C6	C6	D7	D7	E8	E8	E8 ¹
27cm	*	C5	C5	C5	C6	D7	D7	E8	E8	E8	E9 ¹
28cm	*	C5	C5	C5	C6	D7	D7	E8	E8	E9	E9 ¹
29cm	*	*	D6	D6	D7	D7	E8	E8	E9	E9	E9 ¹
30cm	*	*	D6	D6	D7	D7	E8	E8	E9	E9	E9 ¹

* Please contact your Inside Sales Representative for recommendation

¹ For Ibex XD only

Centri® MCV Foot Shell

DF-XXXX13033

	22cm	23cm	24cm	25cm	26cm	27cm	28cm	29cm	30cm
Left	4522	4523	4524	4525	4526	4527	4528	4529	4530
Right	4622	4623	4624	4625	4626	4627	4628	4629	4630

Note: Please allow for a 3 week lead time with custom colored foot shell orders