

» Sport Boot (Model KPB10)

Children's Orthotic Boot

Supplier:

KeepingPace, Inc.

Features & Benefits:

- Inner sole has three varying thicknesses, allowing for easier internal adjustments and customizable depths
- Lightweight
- Ankle support
- Ankle to toe lacing provides a wider opening for easier foot insertion with a custom orthosis or A.F.O.
- Distinctive sole with wider tread design increases direct surface contact and promotes medial/lateral stability
- Wider, deeper heel and toe box
- High abrasion resistant rubber toe cap and reinforced triple stitching for increased durability

Specifications:

- Model: *Black sole, toe and piping*
- Types: *Four (Toddler, Youth, Adolescent and Adult)*
- Sizes: *Toddler= Sizes 7–9*
Youth= Sizes 10–3
Adolescent= Sizes 4–7
- Materials: *High-abrasion rubber, ruva, nylon mesh, leather and non-woven lining materials*
- Do not machine wash or dry—hand wash with a soft cloth dampened with mild soap and cold water

Indications:

- Leg length discrepancies
- Children who wear custom orthoses and A.F.O.s

Call Points:

- Orthotist
- Pedorthist
- Prosthetist

Warranty:

- 30 days from purchase, against manufacturers defects



Sport Boot (Model KPB10) » KeepingPace, Inc.

Item No.	Description	Size (Age)	Size (Foot)*	Color
KEP-KPB10T-7	Lace-up orthotic boot	Toddler	7	Black
KEP-KPB10T-8	Lace-up orthotic boot	Toddler	8	Black
KEP-KPB10T-9	Lace-up orthotic boot	Toddler	9	Black
KEP-KPB10Y-10	Lace-up orthotic boot	Youth	10	Black
KEP-KPB10Y-11	Lace-up orthotic boot	Youth	11	Black
KEP-KPB10Y-12	Lace-up orthotic boot	Youth	12	Black
KEP-KPB10Y-13	Lace-up orthotic boot	Youth	13	Black
KEP-KPB10Y-1	Lace-up orthotic boot	Youth	1	Black
KEP-KPB10Y-2	Lace-up orthotic boot	Youth	2	Black
KEP-KPB10Y-3	Lace-up orthotic boot	Youth	3	Black
KEP-KPB10A-4	Lace-up orthotic boot	Adolescent	4	Black
KEP-KPB10A-5	Lace-up orthotic boot	Adolescent	5	Black
KEP-KPB10A-6	Lace-up orthotic boot	Adolescent	6	Black
KEP-KPB10A-7	Lace-up orthotic boot	Adolescent	7	Black

* Approximate sizes only. For proper fitting, consult with patient or measure their shoe size prior to ordering.