

» HealPRO™

Contracture Boot

Supplier:

Orthomerica Products, Inc.

Features & Benefits:

- Soft CoolFoam™ liner with finished edges with wrap around proximal trim line to reduce calf pressure
- Dual foot and calf Velcro® closures with easy opening tabs to facilitate independent donning and doffing
- Heel offset reduces heel pressure in supine position
- Internal or external rotation control with 90° stop
- Rubber threaded sole to enhance safety
- Toe lift facilitates clearance in swing phase during ambulation
- Extended footplate protects toes from bed linens and prevents contractures
- Universal — fits left or right foot
- Liner is easy to remove and re-apply for cleaning

Specifications:

- Sizes: *Small, Medium, Large and X-Large*
- Straps: *Two*
- Color: *Black with gray lining*
- Extension straps included for edema or bariatric patients
- Optional abductor bar kit is available (sold separately)

Indications:

- Prevention of heel ulcers
- Unloading of pre-existing ulcer on posterior aspect of the calcaneus
- Surgical positioning of the lower limb
- Prevention of lower limb internal or external rotation while patient is supine
- Protection of fragile and compromised skin of the foot and ankle
- Foot drop without contractures
- Patients at risk to lose ankle range of motion secondary to recumbent positioning

Contraindications:

- Fixed ankle contractures
- Excessive swelling

Call Points:

- Orthotist
- Physiotherapist
- Occupational therapist
- General practitioner
- Orthopaedic surgeon
- Rheumatologist

Warranty:

- Against manufacturer defects



Internal or external anti-rotation device with a 90° stop



Adjustable Abductor Bar links to HealPRO when abduction control is required

HealPRO™ » Orthomerica Products, Inc.

Item No.	Description	Size	Max. Calf Circumference*
OR-3530	Foot and ankle positioning recumbent orthosis	Small	16"
OR-3531	Foot and ankle positioning recumbent orthosis	Medium	16"
OR-3532	Foot and ankle positioning recumbent orthosis	Large	18"
OR-3533	Foot and ankle positioning recumbent orthosis	X-Large	20"

Abductor Bar Kit

OR-3529.02	Abductor bar kit for HealPRO™ contracture boot	—	—
------------	--	---	---

HealPRO™ » Sizing Information

Size	Foot Length	Shoe Size (European)	Shoe Size (Women's)	Shoe Size (Men's)
Small	8"– 9"	35–37	3–6	3–5
Medium	9"– 10"	37–41	6–9	5–8
Large	10"– 11"	41–44	9–12	8–11
X-Large	11"– 12"	44–48.5	12–15	11–14