

## » Iceross® Stabilo® Seal-In®

### Liner

#### Supplier:

Ossur Americas

#### Features & Benefits:

- Incorporates a Hypobaric Sealing Membrane (HSM), providing an airtight seal
- Proprietary matrix supports the HSM during donning by preventing elongation of the distal tissue
- Provides the stable and intimate fit required for an active daily routine
- Includes Active Skin Care ingredients to soothe the skin, moisturize, and promote skin health
- Increases range of motion by removing the restriction of a knee sleeve
- Conforms to internal socket shape
- Supplex® cover has non-ravel stitching and superb abrasion tolerance

#### Specifications:

- Type: *Seal-In® Liner*
- Style: *Standard*
- Material: *DermoSil® silicone with Active Skin Care*
- Thickness: *3mm*
- Matrix length: *10cm*
- Sizes: *11 (see chart)*
- When fitting Seal-In liner, a one-way auto expulsion valve with push-button suction release must be used

#### Indications:

- Transtibial amputees
- Moderately active users who require soft tissue control
- Lower limb amputees preferring a suction liner

#### Call Points:

- Prosthetist
- Prosthetic technician

#### Warranty:

- Six months limited, against manufacturers defects



**Iceross® Stabilo® Seal-In® » Ossur Americas**

Item No.	Description	Style	Thickness	Matrix	Size
OS-I-661318	Transfemoral and transtibial Seal-In liner	Standard	3mm	10cm	18cm
OS-I-661320	Transfemoral and transtibial Seal-In liner	Standard	3mm	10cm	20cm
OS-I-661322	Transfemoral and transtibial Seal-In liner	Standard	3mm	10cm	22cm
OS-I-661323	Transfemoral and transtibial Seal-In liner	Standard	3mm	10cm	23.5cm
OS-I-661325	Transfemoral and transtibial Seal-In liner	Standard	3mm	10cm	25cm
OS-I-661326	Transfemoral and transtibial Seal-In liner	Standard	3mm	10cm	26.5cm
OS-I-661328	Transfemoral and transtibial Seal-In liner	Standard	3mm	10cm	28cm
OS-I-661330	Transfemoral and transtibial Seal-In liner	Standard	3mm	10cm	30cm
OS-I-661332	Transfemoral and transtibial Seal-In liner	Standard	3mm	10cm	32cm
OS-I-661334	Transfemoral and transtibial Seal-In liner	Standard	3mm	10cm	34cm
OS-I-661336	Transfemoral and transtibial Seal-In liner	Standard	3mm	10cm	36cm

*Note: When measuring for sizing, please be sure to subtract 1.5cm from measured distal circumference*