

» Rebound™ Air Walker

Foam Semi-pneumatic Walker

Supplier:

Ossur Americas

Features & Benefits:

- Integrated air pump and push-release valve enable tool-free inflation and deflation
- Linked medial and lateral air bladders with single inflation point provide balanced compression — shaped to cradle the achilles, holding the heel in place
- Flex Edge™ material helps prevent the formation of pressure points. Material can also be cut to accommodate larger calf sizes.
- EVA material provides greater shock absorption at the heel
- Full rocker bottom enables patient to have a natural, stable gait
- Trimmable toe bumper provides additional protection
- Ventilated shell provides breathability for increased comfort
- Moisture wicking liner helps keep skin dry and comfortable
- Optional Walker heel wedge device can be used for patients with additional achilles tendon injuries

Specifications:

- Sizes: *X-Small, Small, Medium, Large and X-Large*
- Straps: *Three*
- Color: *Gray*
- Replacement liners are available
- Walker heel wedge is sold separately

Indications:

- Stable fracture of foot and/or ankle
- Severe ankle sprain (Grade II or III)
- Foot and/or ankle trauma
- Post-operative and rehabilitation use

Contraindications:

- Diabetic patients

Call Points:

- Physician
- Orthotist
- Casting technician
- Pedorthist

Warranty:

- Against manufacturer defects



"Images courtesy of Ossur Americas, www.ossur.com"

Rebound™ Air Walker » Ossur Americas

Item No.	Description	Size	Size (Mens)*	Size (Womens)*	Height	Foot Width	Calf Circ.
JM-B-242900001	Foam pneumatic walker	X-Small	2-4	3½-5½	30cm (11 ¹³ / ₁₆ "	10cm (4")	38cm (15")
JM-B-242900002	Foam pneumatic walker	Small	4-7	5½-8½	31.5cm (12 ⁵ / ₃₂ "	12cm (4½")	44cm (17")
JM-B-242900003	Foam pneumatic walker	Medium	7-10½	8½-12	35cm (13¾")	11.4cm (4¾")	51cm (20")
JM-B-242900004	Foam pneumatic walker	Large	10½-12½	12-14	35cm (13¾")	12.7cm (5")	56cm (22")
JM-B-242900005	Foam pneumatic walker	X-Large	12½-15	14-16½	40cm (15¾")	14cm (5½")	58cm (23")

Heel Wedge Insert

JM-B-242902501	Walker heel wedge	Universal			—	—	—
----------------	-------------------	-----------	--	--	---	---	---

Hot and Cold Therapy

JM-WKGEL100	Hot/cold therapy insert	Universal			—	—	—
-------------	-------------------------	-----------	--	--	---	---	---

* Approximate sizes only. For proper fitting, consult with patient or measure their shoe size prior to ordering.

Rebound™ Air Walker » Replacement Liners

Item No.	Description	Size
JM-B-418000101	Rebound™ air walker soft good	X-Small
JM-B-418000102	Rebound™ air walker soft good	Small
JM-B-418000103	Rebound™ air walker soft good	Medium
JM-B-418000104	Rebound™ air walker soft good	Large
JM-B-418000105	Rebound™ air walker soft good	X-Large