

» 600 Classic / MX600 Classic Max

Orthopedic Shoe

Supplier:

Pedors®

Features & Benefits:

- EVA material accomodates a sensitive diabetic foot
- Removable insert allows for accomodation of diabetic inserts or custom orthotics
- Canvas backing provides durability and breathability
- Toe spring promotes smooth gait and reduces the chance of tripping
- Classic Max comes in extra wide widths, and offers a dual strap that allows a wide opening for easy foot insertion

Specifications:

- Sizes: *Classic*= Mens 1, 2, 3 – 13, 14, 15
Womens 3, 4, 5 – 15, 16, 17
Classic Max= Mens 3, 4, 5 – 14, 15
Womens 3, 4, 5 – 16, 17
- Widths: *Classic*= M, W and XW
Classic Max= XW, XXW and XXXW
- Colors: *Classic*= Black and Beige
Classic Max= Black
(Note: Beige color is not available for all sizes. See chart for more information.)

Indications:

- Arthritic patients
- Diabetic patients
- Patients with custom made orthotics
- Patients who require a soft shoe solution for various foot care pathologies
- Patients requiring extra width

Call Points:

- Orthotist
- Chiroprapist
- Podiatrist
- Pedorthist

Warranty:

- Against manufacturers defects



Pedors Classic



Pedors Classic Max

600 Classic / MX600 Classic Max » Pedors®

Item No.	Description	Color	Width	Men's Sizes	Women's Sizes ¹
600 Classic					
PED-600-SS-M	Touch close EVA orthotic shoe	Black	Medium	1, 2, 3 – 13, 14, 15	3, 4, 5 – 15, 16, 17
PED-601-SS-M	Touch close EVA orthotic shoe	Beige	Medium	3½ – 13	5½ – 15
PED-600-SS-W	Touch close EVA orthotic shoe	Black	Wide	1, 2, 3 – 13, 14, 15	3, 4, 5 – 15, 16, 17
PED-601-SS-W	Touch close EVA orthotic shoe	Beige	Wide	3½ – 13	5½ – 15
PED-600-SS-XW	Touch close EVA orthotic shoe	Black	X-Wide	1, 2, 3 – 13, 14, 15	3, 4, 5 – 15, 16, 17
PED-601-SS-XW	Touch close EVA orthotic shoe	Beige	X-Wide	3½ – 13	5½ – 15

MX600 Classic MAX

PED-600MX-SS-XW	Touch close EVA orthotic shoe	Black	X-Wide	3, 4, 5 – 14, 15	3, 4, 5 – 16, 17
PED-600MX-SS-XXW	Touch close EVA orthotic shoe	Black	XX-Wide	3, 4, 5 – 14, 15	3, 4, 5 – 16, 17
PED-600MX-SS-XXXW	Touch close EVA orthotic shoe	Black	XXX-Wide	3, 4, 5 – 14, 15	3, 4, 5 – 16, 17

SS = Size. Note: When ordering a half size, add a .5 after the full shoe size.

ex: PED - 600 - 8.5 - W
 Style 600 Size 8½ Wide width

Important: Use women's sizes when ordering Pedors® footwear. When ordering for men, subtract 2 sizes from women's shoe size.

ex: Women's 13XW = Men's 11XW

600 Classic / MX600 Classic Max » Width Chart

Styles	Width equivalents				
	M	W	XW	XXW	XXXW
Men / Women	B/C	D/E	EE/3E	4E	6E