

» Royal Flush Cool Max® Stretch Sock

Socks

Supplier:

Royal Knit Inc.

Features & Benefits:

- Features CoolMax® fabric
- Quickly moves moisture away from the skin and the breathability of the fiber keeps patient drier and more comfortable
- Round Toe
- Insulates residual limb against the cold
- Form Fitting
- Most sizes fit Symes and PTB patients
- Stretches in length and width
- Latex free
- Color coded to easily identify thickness

Specifications:

- Types: *Standard and Symes*
- Styles: *Standard, with Hole*, Rubberized Hole* or Sewn Edge Top**
(*Please specify when ordering for availability)
- Material: *CoolMax®*
- Thicknesses: *Three (see chart)*
 - 3 ply = *Yellow*
 - 5 ply = *Green*
 - 6 ply = *Gray*
- Sizes: *Seven (see chart)*
- Lengths: *See chart*
- Wash in cold water, using bleach as needed

Indications:

- Protection and comfort of residual limb

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- Against manufacturers defects



Royal Flush Cool Max® Stretch Sock » Royal Knit Inc.

Item No.	Item No.	Item No.	Size	Size	Length	Flat Width	
3 Ply	5 Ply	6 Ply	Product			Top	Toe
Size A (X-Small)							
RK- RF 3-A/B-6	RK- RF 5-A/B-6	RK- RF 6-A/B-6	A	X-Small	4"-6"	4"-4¾"	3"-3½"
RK- RF 3-A/B-10	RK- RF 5-A/B-10	RK- RF 6-A/B-10	A	X-Small	8"-10"	4"-4¾"	3"-3½"
RK- RF 3-A/B-12	RK- RF 5-A/B-12	RK- RF 6-A/B-12	A	X-Small	10"-12"	4"-4¾"	3"-3½"
RK- RF 3-A/B-14	RK- RF 5-A/B-14	RK- RF 6-A/B-14	A	X-Small	12"-14"	4"-4¾"	3"-3½"
RK- RF 3-A/B-18	RK- RF 5-A/B-18	RK- RF 6-A/B-18	A	X-Small	16"-18"	4"-4¾"	3"-3½"
RK- RF 3-A/B-22	RK- RF 5-A/B-22	RK- RF 6-A/B-22	A	X-Small	20"-22"	4"-4¾"	3"-3½"
Size B (X-Small/Small) Note: Will fit PTB Wide-3							
RK- RF 3-B/0-10	RK- RF 5-B/0-10	RK- RF 6-B/0-10	B,0	X-Small/Small	8"-10"	4¾"-6"	3"-3½"
RK- RF 3-B/0-12	RK- RF 5-B/0-12	RK- RF 6-B/0-12	B,0	X-Small/Small	10"-12"	4¾"-6"	3"-3½"
RK- RF 3-B/0-14	RK- RF 5-B/0-14	RK- RF 6-B/0-14	B,0	X-Small/Small	12"-14"	4¾"-6"	3"-3½"
RK- RF 3-B/0-18	RK- RF 5-B/0-18	RK- RF 6-B/0-18	B,0	X-Small/Small	16"-18"	4¾"-6"	3"-3½"
RK- RF 3-B/0-22	RK- RF 5-B/0-22	RK- RF 6-B/0-22	B,0	X-Small/Small	20"-22"	4¾"-6"	3"-3½"
Size 0/1 (Small)							
RK- RF 3-0/1-10	RK- RF 5-0/1-10	RK- RF 6-0/1-10	0,1	Small	8"-10"	6"-7"	3"-4"
RK- RF 3-0/1-12	RK- RF 5-0/1-12	RK- RF 6-0/1-12	0,1	Small	10"-12"	6"-7"	3"-4"
RK- RF 3-0/1-14	RK- RF 5-0/1-14	RK- RF 6-0/1-14	0,1	Small	12"-14"	6"-7"	3"-4"
RK- RF 3-0/1-18	RK- RF 5-0/1-18	RK- RF 6-0/1-18	0,1	Small	16"-18"	6"-7"	3"-4"
RK- RF 3-0/1-22	RK- RF 5-0/1-22	RK- RF 6-0/1-22	0,1	Small	20"-22"	6"-7"	3"-4"
RK- RF 3-0/1-26	RK- RF 5-0/1-26	RK- RF 6-0/1-26	0,1	Small	24"-26"	6"-7"	3"-4"
Size 0/2 (Small/Medium) Note: Will fit PTB Slim-0,1							
RK- RF 3-0/2-10	RK- RF 5-0/2-10	RK- RF 6-0/2-10	0,2	Small/Medium	8"-10"	6"-7"	4"-5"
RK- RF 3-0/2-12	RK- RF 5-0/2-12	RK- RF 6-0/2-12	0,2	Small/Medium	10"-12"	6"-7"	4"-5"
RK- RF 3-0/2-14	RK- RF 5-0/2-14	RK- RF 6-0/2-14	0,2	Small/Medium	12"-14"	6"-7"	4"-5"
RK- RF 3-0/2-18	RK- RF 5-0/2-18	RK- RF 6-0/2-18	0,2	Small/Medium	16"-18"	6"-7"	4"-5"
RK- RF 3-0/2-22	RK- RF 5-0/2-22	RK- RF 6-0/2-22	0,2	Small/Medium	20"-22"	6"-7"	4"-5"
RK- RF 3-0/2-26	RK- RF 5-0/2-26	RK- RF 6-0/2-26	0,2	Small/Medium	24"-26"	6"-7"	4"-5"
Size 1/2 (Medium) Note: VA Standard Size & will fit PTB Slim - 2							
RK- RF 3-1/2-10	RK- RF 5-1/2-10	RK- RF 6-1/2-10	1,2	Medium	8"-10"	7"-8"	4"-5"
RK- RF 3-1/2-12	RK- RF 5-1/2-12	RK- RF 6-1/2-12	1,2	Medium	10"-12"	7"-8"	4"-5"
RK- RF 3-1/2-14	RK- RF 5-1/2-14	RK- RF 6-1/2-14	1,2	Medium	12"-14"	7"-8"	4"-5"
RK- RF 3-1/2-18	RK- RF 5-1/2-18	RK- RF 6-1/2-18	1,2	Medium	16"-18"	7"-8"	4"-5"
RK- RF 3-1/2-22	RK- RF 5-1/2-22	RK- RF 6-1/2-22	1,2	Medium	20"-22"	7"-8"	4"-5"
RK- RF 3-1/2-26	RK- RF 5-1/2-26	RK- RF 6-1/2-26	1,2	Medium	24"-26"	7"-8"	4"-5"

‡ = **Optional** distal opening. Please use the following to indicate if opening is desired:

ex. RK-ORF5-A/B-14

O – Hole

GO – Rubberized Hole

S – Sewn Top Edge

Royal Flush Cool Max® Stretch Sock » Royal Knit Inc.

Item No.	Item No.	Item No.	Size	Size	Length	Flat Width	
3 Ply	5 Ply	6 Ply	Product			Top	Toe
Size 2/3 (Large) Note: VA Standard Size							
RK- RF 3-2/3-10	RK- RF 5-2/3-10	RK- RF 6-2/3-10	2,3	Large	8"-10"	8"-9"	5"-6"
RK- RF 3-2/3-12	RK- RF 5-2/3-12	RK- RF 6-2/3-12	2,3	Large	10"-12"	8"-9"	5"-6"
RK- RF 3-2/3-14	RK- RF 5-2/3-14	RK- RF 6-2/3-14	2,3	Large	12"-14"	8"-9"	5"-6"
RK- RF 3-2/3-18	RK- RF 5-2/3-18	RK- RF 6-2/3-18	2,3	Large	16"-18"	8"-9"	5"-6"
RK- RF 3-2/3-22	RK- RF 5-2/3-22	RK- RF 6-2/3-22	2,3	Large	20"-22"	8"-9"	5"-6"
RK- RF 3-2/3-26	RK- RF 5-2/3-26	RK- RF 6-2/3-26	2,3	Large	24"-26"	8"-9"	5"-6"
Size 3/4 (X-Large) Note: VA Standard Size							
RK- RF 3-3/4-10	RK- RF 5-3/4-10	RK- RF 6-3/4-10	3,4	X-Large	8"-10"	9"-10"	6"-7"
RK- RF 3-3/4-12	RK- RF 5-3/4-12	RK- RF 6-3/4-12	3,4	X-Large	10"-12"	9"-10"	6"-7"
RK- RF 3-3/4-14	RK- RF 5-3/4-14	RK- RF 6-3/4-14	3,4	X-Large	12"-14"	9"-10"	6"-7"
RK- RF 3-3/4-18	RK- RF 5-3/4-18	RK- RF 6-3/4-18	3,4	X-Large	16"-18"	9"-10"	6"-7"
RK- RF 3-3/4-22	RK- RF 5-3/4-22	RK- RF 6-3/4-22	3,4	X-Large	20"-22"	9"-10"	6"-7"
RK- RF 3-3/4-26	RK- RF 5-3/4-26	RK- RF 6-3/4-26	3,4	X-Large	24"-26"	9"-10"	6"-7"
Size 4/5/6 (XX-Large) Note: Will fit PTB Wide – 4,5,6							
RK- RF 3-0/2-10	RK- RF 5-0/2-10	RK- RF 6-0/2-10	4,5,6	XX-Large	8"-10"	10"-12"	7"-8"
RK- RF 3-0/2-12	RK- RF 5-0/2-12	RK- RF 6-0/2-12	4,5,6	XX-Large	10"-12"	10"-12"	7"-8"
RK- RF 3-0/2-14	RK- RF 5-0/2-14	RK- RF 6-0/2-14	4,5,6	XX-Large	12"-14"	10"-12"	7"-8"
RK- RF 3-0/2-18	RK- RF 5-0/2-18	RK- RF 6-0/2-18	4,5,6	XX-Large	16"-18"	10"-12"	7"-8"
RK- RF 3-0/2-22	RK- RF 5-0/2-22	RK- RF 6-0/2-22	4,5,6	XX-Large	20"-22"	10"-12"	7"-8"
RK- RF 3-0/2-26	RK- RF 5-0/2-26	RK- RF 6-0/2-26	4,5,6	XX-Large	24"-26"	10"-12"	7"-8"

‡ = **Optional** distal opening. Please use the following to indicate if opening is desired:

ex. RK-ORF5-A/B-14

O – Hole

GO – Rubberized Hole

S – Sewn Top Edge

Royal Flush Cool Max® Stretch Sock » Royal Knit Inc.

Item No.	Item No.	Item No.	Size	Size	Length	Flat Width	
3 Ply	5 Ply	6 Ply	Product			Top	Toe
Symes Size 0							
RK- RF 3-0/1-18	RK- RF 5-0/1-18	RK- RF 6-0/1-18	0	Small	16"-18" ^{mm}	6"	3"-3½"
RK- RF 3-0/1-22	RK- RF 5-0/1-22	RK- RF 6-0/1-22	0	Small	20"-22"	6"	3"-3½"
RK- RF 3-0/1-26	RK- RF 5-0/1-26	RK- RF 6-0/1-26	0	Small	24"-26"	6"	3"-3½"
RK- RF 3-0/1-18	RK- RF 5-0/1-18	RK- RF 6-0/1-18	0	Small	16"-18"	6"	3"-4"
RK- RF 3-0/1-22	RK- RF 5-0/1-22	RK- RF 6-0/1-22	0	Small	20"-22"	6"	3"-4"
RK- RF 3-0/1-26	RK- RF 5-0/1-26	RK- RF 6-0/1-26	0	Small	24"-26"	6"	3"-4"
RK- RF 3-0/1-18	RK- RF 5-0/1-18	RK- RF 6-0/1-18	0	Small	16"-18"	6"	3½"-4"
RK- RF 3-0/1-22	RK- RF 5-0/1-22	RK- RF 6-0/1-22	0	Small	20"-22"	6"	3½"-4"
RK- RF 3-0/1-26	RK- RF 5-0/1-26	RK- RF 6-0/1-26	0	Small	24"-26"	6"	3½"-4"
RK- RF 3-0/2-18	RK- RF 5-0/2-18	RK- RF 6-0/2-18	0	Small/Medium	16"-18"	6"	3½"-4½"
RK- RF 3-0/2-22	RK- RF 5-0/2-22	RK- RF 6-0/2-22	0	Small/Medium	20"-22"	6"	3½"-4½"
RK- RF 3-0/2-26	RK- RF 5-0/2-26	RK- RF 6-0/2-26	0	Small/Medium	24"-26"	6"	3½"-4½"
Symes Size 1							
RK- RF 3-0/1-18	RK- RF 5-0/1-18	RK- RF 6-0/1-18	1	Small	16"-18" ^{mm}	7"	3"-3½"
RK- RF 3-0/1-22	RK- RF 5-0/1-22	RK- RF 6-0/1-22	1	Small	20"-22"	7"	3"-3½"
RK- RF 3-0/1-26	RK- RF 5-0/1-26	RK- RF 6-0/1-26	1	Small	24"-26"	7"	3"-3½"
RK- RF 3-0/1-18	RK- RF 5-0/1-18	RK- RF 6-0/1-18	1	Small	16"-18"	7"	3"-4"
RK- RF 3-0/1-22	RK- RF 5-0/1-22	RK- RF 6-0/1-22	1	Small	20"-22"	7"	3"-4"
RK- RF 3-0/1-26	RK- RF 5-0/1-26	RK- RF 6-0/1-26	1	Small	24"-26"	7"	3"-4"
RK- RF 3-0/1-18	RK- RF 5-0/1-18	RK- RF 6-0/1-18	1	Small	16"-18"	7"	3½"-4"
RK- RF 3-0/1-22	RK- RF 5-0/1-22	RK- RF 6-0/1-22	1	Small	20"-22"	7"	3½"-4"
RK- RF 3-0/1-26	RK- RF 5-0/1-26	RK- RF 6-0/1-26	1	Small	24"-26"	7"	3½"-4"
RK- RF 3-1/2-18	RK- RF 5-1/2-18	RK- RF 6-1/2-18	1	Medium	16"-18"	7"	3½"-4½"
RK- RF 3-1/2-22	RK- RF 5-1/2-22	RK- RF 6-1/2-22	1	Medium	20"-22"	7"	3½"-4½"
RK- RF 3-1/2-26	RK- RF 5-1/2-26	RK- RF 6-1/2-26	1	Medium	24"-26"	7"	3½"-4½"

‡ = **Optional** distal opening. Please use the following to indicate if opening is desired:

ex. RK-ORF5-A/B-14

O – Hole

GO – Rubberized Hole

S – Sewn Top Edge