

Synthetic Slipper Casting Sock

Casting Application

Supplier:

STS Company

Features & Benefits:

- Ease of use
- Comfortable and accurate
- Clean — no plaster mess
- Durable — virtually unbreakable once cured
- Reduces casting time for foot orthotics by 50%
- You can now obtain a superior foot mold in less than 2 minutes
- No more waiting until the next day to ship the casts to your lab
- Knitted slipper shaped sock
- Cast can be taken off foot or inside of shoe gear when the baggie is applied to the outside of the slipper sock

Specifications:

- Slipper sock
- Ultra-fast setting resin
- Made from polyester
- Knitted and then impregnated with a water curable polyurethane resin

Indications:

- Custom functional foot orthotic devices
- Custom accommodative foot orthotic devices
- Custom shoe insoles
- Post-op immobilization (applying a double layer)

Call Points:

- Orthotist
- Orthotic technician
- Pedorthist

Warranty:

- Against manufacturer defects



See Materials Catalog 2007/08, #22 (page 25)

Item No.	Description	Size ►	Women	Men	Unit
STS-903-S	Synthetic slipper casting sock	Small	5–6.5	–	10 pairs/Box
STS-903-M	Synthetic slipper casting sock	Medium	7–10	7–9	10 pairs/Box
STS-903-L	Synthetic slipper casting sock	Large	10.5–12	9.5–11.5	10 pairs/Box
STS-903-XL	Synthetic slipper casting sock	X-Large	–	12–15	10 pairs/Box