

## » Advance™ Comfort Laced Footwear

### Orthopedic Shoe

#### Supplier:

Soletech, Inc.

#### Features & Benefits:

- Insta-Flex® sole construction for lightweight EVA cushioning
- High abrasion rubber outsole with high slip resistance
- Full grain leather uppers
- New Dri-Lex® Bamboo Eco moisture wicking / anti-bacterial lining
- Provides medium depth last for a wide range of feet
- Deep toe box
- Removeable insoles to accomodate custom othotics and AFO's

#### Specifications:

- Sizes: *Mens* = 6 – 8½ (M, W and XW)  
9 – 12, 13 and 14 (N, M, W and XW)  
*Womens* = 5 – 6½ (M, W and XW)  
7 – 10, 11 and 12 (N, M, W and XW)
- Widths: *Mens* = N (B/C), M (D), W (E/EE), XW (3E/4E)  
*Womens* = N (A), M (B/C), W (E), XW (3E)
- Depth: *Medium*
- Nylon laces
- Comes with 2 insoles for ease in adjusting the fit

#### Indications:

- Diabetic patients
- Patients with custom made orthotics
- Patients requiring extra depth

#### Call Points:

- Orthotist
- Chiroprapist
- Podiatrist
- Pedorthist

#### Warranty:

- Against manufacturers defects



*Men's black smooth leather shoe*



*Women's light mocha smooth leather shoe*

**Advance™ Comfort Laced Footwear » Soletech, Inc.**

Item No.	Description	Color	Width	Sizes
<b>Men's Footwear</b>				
SOL-E3210-SSN	Men's smooth leather orthotic shoe	Black	Narrow	9 – 12, 13 & 14
SOL-E3210-SSM	Men's smooth leather orthotic shoe	Black	Medium	6 – 12, 13 & 14
SOL-E3210-SSW	Men's smooth leather orthotic shoe	Black	Wide	6 – 12, 13 & 14
SOL-E3210-SSXW	Men's smooth leather orthotic shoe	Black	Extra Wide	6 – 12, 13 & 14
SOL-E3225-SSN	Men's smooth leather orthotic shoe	Rust	Narrow	9 – 12, 13 & 14
SOL-E3225-SSM	Men's smooth leather orthotic shoe	Rust	Medium	6 – 12, 13 & 14
SOL-E3225-SSW	Men's smooth leather orthotic shoe	Rust	Wide	6 – 12, 13 & 14
SOL-E3225-SSXW	Men's smooth leather orthotic shoe	Rust	Extra Wide	6 – 12, 13 & 14

<b>Women's Footwear</b>				
SOL-E3710-SSN	Women's smooth leather orthotic shoe	Black	Narrow	7 – 10, 11 & 12
SOL-E3710-SSM	Women's smooth leather orthotic shoe	Black	Medium	5 – 10, 11 & 12
SOL-E3710-SSW	Women's smooth leather orthotic shoe	Black	Wide	5 – 10, 11 & 12
SOL-E3710-SSXW	Women's smooth leather orthotic shoe	Black	Extra Wide	5 – 10, 11 & 12
SOL-E3736-SSN	Women's smooth leather orthotic shoe	Bone	Narrow	7 – 10, 11 & 12
SOL-E3736-SSM	Women's smooth leather orthotic shoe	Bone	Medium	5 – 10, 11 & 12
SOL-E3736-SSW	Women's smooth leather orthotic shoe	Bone	Wide	5 – 10, 11 & 12
SOL-E3736-SSXW	Women's smooth leather orthotic shoe	Bone	Extra Wide	5 – 10, 11 & 12
SOL-E3744-SSN	Women's smooth leather orthotic shoe	Light Mocha	Narrow	7 – 10, 11 & 12
SOL-E3744-SSM	Women's smooth leather orthotic shoe	Light Mocha	Medium	5 – 10, 11 & 12
SOL-E3744-SSW	Women's smooth leather orthotic shoe	Light Mocha	Wide	5 – 10, 11 & 12
SOL-E3744-SSXW	Women's smooth leather orthotic shoe	Light Mocha	Extra Wide	5 – 10, 11 & 12

SS = Size. Note: When ordering a half size, add a 5 after the full shoe size.

ex: SOL - E3210 - 85W  
Style E-3210 8½ Wide

**Advance™ Footwear » Width Chart**

Styles	Width equivalents			
	N	M	W	XW
Men	B/C	D	E/EE	3E/4E
Women	A	B/C	E	3E