

» Advance™ Comfort Velcro® Footwear

Orthopedic Shoe

Supplier:

Soletech, Inc.

Features & Benefits:

- Insta-Flex® sole construction for lightweight EVA cushioning
- High abrasion rubber outsole with high slip resistance
- Full grain leather uppers
- New Dri-Lex® Bamboo Eco moisture wicking / anti-bacterial lining
- Provides medium depth last for a wide range of feet
- Deep toe box
- Removeable insoles to accomodate custom othotics and AFO's

Specifications:

- Sizes: *Mens* = 6 – 8½ (M, W and XW)
9 – 12, 13 and 14 (N, M, W and XW)
Womens = 5 – 6½ (M, W and XW)
7 – 10, 11 and 12 (N, M, W and XW)
- Widths: *Mens* = N (B/C), M (D), W (E/EE), XW (3E/4E)
Womens = N (A), M (B/C), W (E), XW (3E)
- Depth: *Medium*
- Velcro® straps
- Comes with 2 insoles for ease in adjusting the fit

Indications:

- Diabetic patients
- Patients with custom made orthotics
- Patients requiring extra depth

Call Points:

- Orthotist
- Chiroprapist
- Podiatrist
- Pedorthist

Warranty:

- Against manufacturers defects



Men's rust Velcro® smooth leather shoe



Women's light mocha smooth leather shoe

Advance™ Comfort Velcro® Footwear » Soletech, Inc.

Item No.	Description	Color	Width	Sizes
Men's Footwear				
SOL-E3310-SSN	Men's Velcro® smooth leather orthotic shoe	Black	Narrow	9 – 12, 13 & 14
SOL-E3310-SSM	Men's Velcro® smooth leather orthotic shoe	Black	Medium	6 – 12, 13 & 14
SOL-E3310-SSW	Men's Velcro® smooth leather orthotic shoe	Black	Wide	6 – 12, 13 & 14
SOL-E3310-SSXW	Men's Velcro® smooth leather orthotic shoe	Black	Extra Wide	6 – 12, 13 & 14
SOL-E3325-SSN	Men's Velcro® smooth leather orthotic shoe	Rust	Narrow	9 – 12, 13 & 14
SOL-E3325-SSM	Men's Velcro® smooth leather orthotic shoe	Rust	Medium	6 – 12, 13 & 14
SOL-E3325-SSW	Men's Velcro® smooth leather orthotic shoe	Rust	Wide	6 – 12, 13 & 14
SOL-E3325-SSXW	Men's Velcro® smooth leather orthotic shoe	Rust	Extra Wide	6 – 12, 13 & 14

Women's Footwear

SOL-E3810-SSN	Women's Velcro® smooth leather orthotic shoe	Black	Narrow	7 – 10, 11 & 12
SOL-E3810-SSM	Women's Velcro® smooth leather orthotic shoe	Black	Medium	5 – 10, 11 & 12
SOL-E3810-SSW	Women's Velcro® smooth leather orthotic shoe	Black	Wide	5 – 10, 11 & 12
SOL-E3810-SSXW	Women's Velcro® smooth leather orthotic shoe	Black	Extra Wide	5 – 10, 11 & 12
SOL-E3836-SSN	Women's Velcro® smooth leather orthotic shoe	Bone	Narrow	7 – 10, 11 & 12
SOL-E3836-SSM	Women's Velcro® smooth leather orthotic shoe	Bone	Medium	5 – 10, 11 & 12
SOL-E3836-SSW	Women's Velcro® smooth leather orthotic shoe	Bone	Wide	5 – 10, 11 & 12
SOL-E3836-SSXW	Women's Velcro® smooth leather orthotic shoe	Bone	Extra Wide	5 – 10, 11 & 12
SOL-E3844-SSN	Women's Velcro® smooth leather orthotic shoe	Light Mocha	Narrow	7 – 10, 11 & 12
SOL-E3844-SSM	Women's Velcro® smooth leather orthotic shoe	Light Mocha	Medium	5 – 10, 11 & 12
SOL-E3844-SSW	Women's Velcro® smooth leather orthotic shoe	Light Mocha	Wide	5 – 10, 11 & 12
SOL-E3844-SSXW	Women's Velcro® smooth leather orthotic shoe	Light Mocha	Extra Wide	5 – 10, 11 & 12

SS = Size. Note: When ordering a half size, add a 5 after the full shoe size.

ex: SOL - E3310 - 85W
 Style E-3310 8½ Wide

Advance™ Footwear » Width Chart

Styles	Width equivalents			
	N	M	W	XW
Men	B/C	D	E/EE	3E/4E
Women	A	B/C	E	3E