

**Epi-med® Plus**  
**Epicondylitis Clasp**

**Vendor:**

Thämert GmbH & Co.

**Features & Benefits:**

- Clasp compresses the forearm musculature at its largest circumference below the elbow
- Compression achieved through two pads, not a circular bandage, so will not impair bloodflow
- Lightweight
- Removable felt insert incorporates silicone pads
- Strap end with a white arrow allows for visual check of tension control
- Blue silicone inserts minimize slippage

**Specifications:**

- Sized by measuring the forearm approximately 4–5 cm (1½"–2") below the elbow
- Available in five sizes (see table)
- Universal fit: left or right

**Indications:**

- Tendopatias in the elbow region, caused by muscular overuse after work, sports or trauma
- Epicondylitis humeri lateralis (tennis elbow)
- Epicondylitis humeri medialis (golfer's elbow)

**Contraindications:**

- Osteochondrosis

**Call Points:**

- Orthotist
- Physiotherapist
- Sports medicine therapist
- Occupational therapist
- Athletic trainer

**Warranty:**

- Against manufacturer defects



1767

**Epi-med® Plus – Thämert GmbH & Co.**

Item No.	Description	Size	Arm Circumference	Color
TH-48/007-XS	Epi-med® Plus Clasp	X-Small	8½"–9" (22–23cm)	White w/blue trim
TH-48/007-SM	Epi-med® Plus Clasp	Small	9¼"–10" (24–25cm)	White w/blue trim
TH-48/007-MD	Epi-med® Plus Clasp	Medium	10¼"–10¾" (26–27cm)	White w/blue trim
TH-48/007-LG	Epi-med® Plus Clasp	Large	11"–11¾" (28–30cm)	White w/blue trim
TH-48/007-XL	Epi-med® Plus Clasp	X-Large	12"–13½" (31–34cm)	White w/blue trim
TH-48/017-XS	Replacement Felt Insert w/Silicone Pad	X-Small	8½"–9" (22–23cm)	White w/blue silicone
TH-48/017-SM	Replacement Felt Insert w/Silicone Pad	Small	9¼"–10" (24–25cm)	White w/blue silicone
TH-48/017-MD	Replacement Felt Insert w/Silicone Pad	Medium	10¼"–10¾" (26–27cm)	White w/blue silicone
TH-48/017-LG	Replacement Felt Insert w/Silicone Pad	Large	11"–11¾" (28–30cm)	White w/blue silicone
TH-48/017-XL	Replacement Felt Insert w/Silicone Pad	X-Large	12"–13½" (31–34cm)	White w/blue silicone