

Thermocork

Vendor:
OrtoPed

Applications:
• For complex and lasting connections

Patient Benefits:
• Durable

Practitioner Benefits:
• Can be formed directly on positive model
• Easily bonds to other materials

Features:

- Thermoformable
- Adheres easily when heated
- Resistant to compression — will maintain its shape even under severe stress, does not bottom out
- Strong but flexible
- Easy to grind and trim
- Moldable and vacuum formable
- Localized pressure points may be compressed with use of heat gun

Specifications:

- Combined natural cork particles
- Durometer Shore A: 55 ±5
- Heating temperature: 170°F – 200°F

Indications:

- Total contact orthotics
- Patients where weight is a factor
- Modifications
- Posting

Call Points:

- Orthotist
- Orthotic technician
- Pedorthist



Thermocork

Item No.	Sheet Size
APX-8013112	1/8" x 36" x 44"
APX-8013113	3/16" x 36" x 44"
APX-8013114	1/4" x 36" x 44"
APX-8013115	3/8" x 36" x 44"