

» Seattle Carbon LightFoot®

Prosthetic Foot

Supplier:

Trulife Group Ltd.

Features & Benefits:

- Unique split composite keel
- Increased dynamic response and ML compliance
- Increased stability, particularly on uneven terrain and during high activity
- Available in soft and firm foam durometers
- Three different keels available to accommodate patient's weight and activity level

Specifications:

- Types: *Standard and Ladies*
- Sizes: *Standard= 22–29 cm*
Ladies= 22–28 cm
- Weight limit: *125kg (275lbs)*
- Weight: *373g (13.2oz) (Based on size 25cm standard foot)*
- Keel stiffnesses: *1–3 (see chart)*
- Heel height: *10mm (3/8")*
- Clearance: *38mm (1 1/2")*
- Foot shell: *Caucasian, Light Brown or Dark Brown*
- Activity level: *Low to Moderate (K1–K2)*
(Note: Patient weight determines optimal activity level)
- Includes 10M x 55 metric bolt kit (SF-SFB540)

Indications:

- K1–K2 transfemoral and transtibial amputees

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- 24 months on keel from date of installation, against manufacturers defects
- 12 months on the foot shell, against manufacturers defects



Seattle Carbon LightFoot® » Trulife Group Ltd.

Ordering example: SF-SCF185-24L-N3

Item No.	Description	Side	Color	Clearance	Heel Height	Weight ¹
Standard						
SF-SCF185-SSL-NK	Prosthetic foot	Left	Caucasian	38mm (1½")	10mm (¾")	373g (13.2oz)
SF-SCF185-SSR-NK	Prosthetic foot	Right	Caucasian	38mm (1½")	10mm (¾")	373g (13.2oz)
SF-SCF186-SSL-NK	Prosthetic foot	Left	Light Brown	38mm (1½")	10mm (¾")	373g (13.2oz)
SF-SCF186-SSR-NK	Prosthetic foot	Right	Light Brown	38mm (1½")	10mm (¾")	373g (13.2oz)
SF-SCF188-SSL-NK	Prosthetic foot	Left	Dark Brown	38mm (1½")	10mm (¾")	373g (13.2oz)
SF-SCF188-SSR-NK	Prosthetic foot	Right	Dark Brown	38mm (1½")	10mm (¾")	373g (13.2oz)

Ladies

SF-SCF195-SSL-NK	Ladies prosthetic foot	Left	Caucasian	38mm (1½")	10mm (¾")	373g (13.2oz)
SF-SCF195-SSR-NK	Ladies prosthetic foot	Right	Caucasian	38mm (1½")	10mm (¾")	373g (13.2oz)
SF-SCF196-SSL-NK	Ladies prosthetic foot	Left	Light Brown	38mm (1½")	10mm (¾")	373g (13.2oz)
SF-SCF196-SSR-NK	Ladies prosthetic foot	Right	Light Brown	38mm (1½")	10mm (¾")	373g (13.2oz)
SF-SCF198-SSL-NK	Ladies prosthetic foot	Left	Dark Brown	38mm (1½")	10mm (¾")	373g (13.2oz)
SF-SCF198-SSR-NK	Ladies prosthetic foot	Right	Dark Brown	38mm (1½")	10mm (¾")	373g (13.2oz)

Foot Bolt

SF-SFB540 Foot bolt (10M x 55mm)

SS= Size (22, 24, 25, 26, 27, 28, 29², 30²)

K= Keel (1-3)

¹ Based on a 25cm standard foot

² Sizes available for standard type only

Seattle Carbon LightFoot® » Keel Selection (K)

Keel	Low (Walking with Aid)	Medium-Low (Limited Walking)	Medium (Walking)	Medium-High (Jogging, Light Sports)	High (Running, Basketball)
1	< 68 kg (< 150lbs)	< 54kg (< 120lbs)	< 45kg (< 100lbs)	< 38kg (< 85lbs)	< 34kg (< 75lbs)
2	68–113 kg (150–250 lbs)	54–95 kg (120–210 lbs)	45–79 kg (100–175 lbs)	38–65 kg (85–145 lbs)	34–56 kg (75–125 lbs)
3	114–125 kg (251–275 lbs)	96–125 kg (211–275 lbs)	80–122 kg (176–270 lbs)	66–106 kg (146–235 lbs)	57–90 kg (126–220 lbs)

Note: If the patient has a long transtibial amputation, carries heavy loads or will reach a higher activity level within a year, then select a stiffer keel.

Important: Selecting a keel that is less stiff than the above recommendations will void the warranty.