

» Seattle Lightfoot² with Titanium Pyramid

Prosthetic Foot

Supplier:

Trulife Group Ltd.

Features & Benefits:

- Integrated titanium pyramid
- Unique split keel
- Increased dynamic response and ML compliance
- Increased stability, particularly on uneven terrain

Specifications:

- Sizes: 22–30 cm*
(* 30cm size only available for 19mm heel height models)
- Weight limit: 136kg (300lbs)
- Weight: 480g (16.9oz) (Based on size 25cm foot)
- Keel stiffnesses: 5–9 (see chart)
- Heel height: 10mm ($\frac{3}{8}$ ") and 19mm ($\frac{3}{4}$ ")
- Clearance: 90mm ($3\frac{1}{2}$ ")
- Foot shell: *Caucasian, Light Brown or Dark Brown*
- Activity level: *Moderate (K2)*
(Note: Patient weight determines optimal activity level)
- Includes 10M x 80 metric bolt kit (SF-SFB542), 25mm foot spacer (SF-SFS232) and foot pyramid (SF-SCA240)

Indications:

- K2 Transfemoral and transtibial amputees

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- 24 months on keel from date of installation, against manufacturers defects
- 12 months on the foot shell, against manufacturers defects



Seattle Lightfoot² with Titanium Pyramid » Trulife Group Ltd.

Ordering example: SF-SLF190-24L-H6

Item No.	Description	Side	Color	Clearance	Heel Height	Weight ¹
SF-SLF195-SSL-HK	Prosthetic foot	Left	Caucasian	90mm (3½")	10mm (¾")	480g (16.9oz)
SF-SLF195-SSR-HK	Prosthetic foot	Right	Caucasian	90mm (3½")	10mm (¾")	480g (16.9oz)
SF-SLF196-SSL-HK	Prosthetic foot	Left	Light Brown	90mm (3½")	10mm (¾")	480g (16.9oz)
SF-SLF196-SSR-HK	Prosthetic foot	Right	Light Brown	90mm (3½")	10mm (¾")	480g (16.9oz)
SF-SLF198-SSL-HK	Prosthetic foot	Left	Dark Brown	90mm (3½")	10mm (¾")	480g (16.9oz)
SF-SLF198-SSR-HK	Prosthetic foot	Right	Dark Brown	90mm (3½")	10mm (¾")	480g (16.9oz)
SF-SLF190-SSL-GK	Prosthetic foot	Left	Caucasian	90mm (3½")	19mm (¾")	480g (16.9oz)
SF-SLF190-SSR-GK	Prosthetic foot	Right	Caucasian	90mm (3½")	19mm (¾")	480g (16.9oz)
SF-SLF191-SSL-GK	Prosthetic foot	Left	Light Brown	90mm (3½")	19mm (¾")	480g (16.9oz)
SF-SLF191-SSR-GK	Prosthetic foot	Right	Light Brown	90mm (3½")	19mm (¾")	480g (16.9oz)
SF-SLF193-SSL-GK	Prosthetic foot	Left	Dark Brown	90mm (3½")	19mm (¾")	480g (16.9oz)
SF-SLF193-SSR-GK	Prosthetic foot	Right	Dark Brown	90mm (3½")	19mm (¾")	480g (16.9oz)

SS= Size (22, 24, 25, 26, 27, 28, 29, 30²)

K= Keel (5–9)

¹ Based on a 25cm foot

² Size available only for 19mm heel height

Seattle Lightfoot² with Titanium Pyramid » Keel Selection (K)

Keel	Low (Walking with Aid)	Medium-Low (Limited Walking)	Medium (Walking)	Medium-High (Jogging, Light Sports)	High (Running, Basketball)
5	<82kg (<181lbs)	<59kg (<131lbs)	<45kg (<101lbs)	<36kg (<81lbs)	<29kg (<66lbs)
6	82–122 kg (181–270 lbs)	59–74 kg (131–165 lbs)	45–67 kg (101–150 lbs)	36–53 kg (81–120 lbs)	29–44 kg (66–100 lbs)
7	123–136 kg (271–300 lbs)	75–85 kg (166–190 lbs)	68–90 kg (151–200 lbs)	54–72 kg (121–160 lbs)	45–60 kg (101–135 lbs)
8	—	86–115 kg (191–255 lbs)	91–113 kg (201–250 lbs)	73–90 kg (161–200 lbs)	61–76 kg (136–169 lbs)
9	—	116–136 kg (256–300 lbs)	114–136 kg (251–300 lbs)	91–106 kg (201–235 lbs)	77–90 kg (170–200 lbs)

Note: If the patient has a long transtibial amputation, carries heavy loads or will reach a higher activity level within a year, then select a stiffer keel.

Important: Selecting a keel that is less stiff than the above recommendations will void the warranty.