

» Seattle Turas

Prosthetic Foot

Supplier:

Trulife Group, Ltd.

Features & Benefits:

- Integrated multi-axial composite ankle
- Dynamic split keel design
- Enhanced plantar flexion at heel strike provides early foot stability

Specifications:

- Sizes: 22–29 cm
- Weight limit: 113kg (250lbs)
- Weight: 560g (19.8oz) (based on size 25cm foot, without foot shell)
- Keel stiffnesses: 5–8 (see chart)
- Heel height: 10mm ($\frac{3}{8}$ ")
- Clearance: 85mm ($3\frac{3}{8}$ ")
- Foot shell: Split toe foot shell in Caucasian, Light Brown or Dark Brown
- Activity level: High (K3)
(Note: Patient weight determines optimal activity level)
- Foot shell included

Indications:

- Active K3 transfemoral and transtibial amputees

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- 24 months on keel from date of installation, against manufacturers defects
- 6 months on foot shell, against manufacturers defects



Seattle Turas foot without foot shell

Seattle Turas » Trulife Group, Ltd.

Ordering example: SF-SLP300-24-R-J6

Item No.	Description	Side	Shell Color	Clearance	Heel Height	Weight ¹
SF-SLP300-SS-L-JKX	High-impact prosthetic foot	Left	Caucasian	85mm (3 ³ / ₈ "	10mm (3 ³ / ₈ "	560g (19.8oz)
SF-SLP300-SS-R-JKX	High-impact prosthetic foot	Right	Caucasian	85mm (3 ³ / ₈ "	10mm (3 ³ / ₈ "	560g (19.8oz)
SF-SLP301-SS-L-JKX	High-impact prosthetic foot	Left	Light Brown	85mm (3 ³ / ₈ "	10mm (3 ³ / ₈ "	560g (19.8oz)
SF-SLP301-SS-R-JKX	High-impact prosthetic foot	Right	Light Brown	85mm (3 ³ / ₈ "	10mm (3 ³ / ₈ "	560g (19.8oz)
SF-SLP303-SS-L-JKX	High-impact prosthetic foot	Left	Dark Brown	85mm (3 ³ / ₈ "	10mm (3 ³ / ₈ "	560g (19.8oz)
SF-SLP303-SS-R-JKX	High-impact prosthetic foot	Right	Dark Brown	85mm (3 ³ / ₈ "	10mm (3 ³ / ₈ "	560g (19.8oz)

SS= Size (22, 23, 24, 25, 26, 27, 28, 29)

K= Keel (5, 6, 7, 8)

X= Size code (2, 3, 4, 5, 6, 7, 8, 9). See below for more information.

¹ Based on a 25cm foot, without foot shell

Seattle Turas » Keel Stiffness (K)

Keel	Low (Walking with aid)	Medium-Low (Limited Walking)	Medium (Walking on Uneven Surfaces)	Medium-High (Light Sports)	High (Running, Basketball)
5	< 81kg < 180lbs	< 59kg < 130lbs	< 45kg < 100lbs	< 36kg < 80lbs	< 29kg < 65lbs
6	81–113 kg (180–250 lbs)	59–86 kg (130–190 lbs)	45–68 kg (100–150 lbs)	36–54 kg (80–120 lbs)	29–45 kg (65–100 lbs)
7	—	87–113 kg (191–250 lbs)	69–90 kg (151–200 lbs)	55–72 kg (121–160 lbs)	46–61 kg (101–135 lbs)
8	—	—	91–113 kg (201–250 lbs)	73–90 kg (161–200 lbs)	62–74 kg (136–165 lbs)

Note: If the patient has a long transtibial amputation, carries heavy loads or will reach a higher activity level within a year, then select a stiffer keel.

Important: Selecting a keel that is less stiff than the above recommendations will void the warranty.

Seattle Turas » Size Code (X)

Size	22cm	23cm	24cm	25cm	26cm	27cm	28cm	29cm
Code	2	3	4	5	6	7	8	9

Seattle Turas » Parts

Item No.	Description	Shell color	Side
Foot Shell			
SF-SFC300-SS-L	Foot shell for Seattle Turas	Caucasian	Left
SF-SFC300-SS-R	Foot shell for Seattle Turas	Caucasian	Right
SF-SFC301-SS-L	Foot shell for Seattle Turas	Light Brown	Left
SF-SFC301-SS-R	Foot shell for Seattle Turas	Light Brown	Right
SF-SFC303-SS-L	Foot shell for Seattle Turas	Dark Brown	Left
SF-SFC303-SS-R	Foot shell for Seattle Turas	Dark Brown	Right
Tendon Replacement Kit			
SF-SLP300-KIT	Tendon replacement kit for Seattle Turas (includes: Hex nut, Wire and Tendon)		