

» Carbon Copy® High Performance

Prosthetic Foot

Supplier:

WillowWood™

Features & Benefits:

- Designed for moderate to highly active amputees
- Unique design responds to both low and high activity levels
- Keel/plate attachment allows energy storage at mid-stance
- Molded titanium pyramid reduces the foot's overall weight
- Available in three heel heights

Specifications:

- Sizes: *Mens= 23–31 cm*
Womens= 22–27 cm
- Weight limit: *113kg (250lbs)*
- Weight: *Mens= 635g (22.3oz)*
Womens= 585g (20.7oz)
(Based on 26cm size foot)
- Toe spring stiffnesses: *1–5 (see chart)*
- Heel spring stiffnesses: *Soft, Medium or Firm*
- Heel heights: *10mm (3/8"), 19mm (3/4)"* or 30mm (1 1/8)"***
*(*Mens foot only)*
*(**Womens foot only)*
- Clearance: *89mm (3 1/2")*
- Foot shell: *Unisex, Lifelike in Buff, Tan or Medium Brown*
- Foot shell height: *76mm (3")*
- Activity level: *High to Very High (K3-K4)*
(Note: Patient weight determines optimal activity level)
- All heel inserts are included

Indications:

- Transtibial and transfemoral amputees
- Active amputees

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- One year and three weeks from purchase, against manufacturers defects



Carbon Copy® High Performance » WillowWood™

Ordering example: OW-CHPM-24RB-3L

Item No.	Description	Gender	Side	Color	Heel Height	Weight ³
OW-CHPM-SSLB-KL	High activity prosthetic foot	Mens	Left	Buff	10mm (3/8")	635g (22.4oz)
OW-CHPM-SSRB-KL	High activity prosthetic foot	Mens	Right	Buff	10mm (3/8")	635g (22.4oz)
OW-CHPW-SSLB-KL	High activity prosthetic foot	Womens	Left	Buff	10mm (3/8")	585g (20.6oz)
OW-CHPW-SSRB-KL	High activity prosthetic foot	Womens	Right	Buff	10mm (3/8")	585g (20.6oz)
OW-CHPM-SSLM-KL	High activity prosthetic foot	Mens	Left	Medium brown	10mm (3/8")	635g (22.4oz)
OW-CHPM-SSRM-KL	High activity prosthetic foot	Mens	Right	Medium brown	10mm (3/8")	635g (22.4oz)
OW-CHPW-SSLM-KL	High activity prosthetic foot	Womens	Left	Medium brown	10mm (3/8")	585g (20.6oz)
OW-CHPW-SSRM-KL	High activity prosthetic foot	Womens	Right	Medium brown	10mm (3/8")	585g (20.6oz)
OW-CHPM-SSLT-KL	High activity prosthetic foot	Mens	Left	Tan	10mm (3/8")	635g (22.4oz)
OW-CHPM-SSRT-KL	High activity prosthetic foot	Mens	Right	Tan	10mm (3/8")	635g (22.4oz)
OW-CHPW-SSLT-KL	High activity prosthetic foot	Womens	Left	Tan	10mm (3/8")	585g (20.6oz)
OW-CHPW-SSRT-KL	High activity prosthetic foot	Womens	Right	Tan	10mm (3/8")	585g (20.6oz)
OW-CHPM-SSLB-KM	High activity prosthetic foot	Mens	Left	Buff	19mm (3/4")	635g (22.4oz)
OW-CHPM-SSRB-KM	High activity prosthetic foot	Mens	Right	Buff	19mm (3/4")	635g (22.4oz)
OW-CHPM-SSLM-KM	High activity prosthetic foot	Mens	Left	Medium brown	19mm (3/4")	635g (22.4oz)
OW-CHPM-SSRM-KM	High activity prosthetic foot	Mens	Right	Medium brown	19mm (3/4")	633g (22.3oz)
OW-CHPM-SSLT-KM	High activity prosthetic foot	Mens	Left	Tan	19mm (3/4")	635g (22.4oz)
OW-CHPM-SSRT-KM	High activity prosthetic foot	Mens	Right	Tan	19mm (3/4")	635g (22.3oz)
OW-CHPW-SSLB-KH	High activity prosthetic foot	Womens	Left	Buff	30mm (1 1/8")	585g (20.6oz)
OW-CHPW-SSRB-KH	High activity prosthetic foot	Womens	Right	Buff	30mm (1 1/8")	585g (20.6oz)
OW-CHPW-SSLM-KH	High activity prosthetic foot	Womens	Left	Medium brown	30mm (1 1/8")	585g (20.6oz)
OW-CHPW-SSRM-KH	High activity prosthetic foot	Womens	Right	Medium brown	30mm (1 1/8")	585g (20.6oz)
OW-CHPW-SSLT-KH	High activity prosthetic foot	Womens	Left	Tan	30mm (1 1/8")	585g (20.6oz)
OW-CHPW-SSRT-KH	High activity prosthetic foot	Womens	Right	Tan	30mm (1 1/8")	585g (20.6oz)

SS = Size (22², 23, 24, 25, 26, 27, 28¹, 29¹, 30¹, 31¹)

K = Toe Spring (1–5)

¹ For mens styles only

² For womens styles only

³ Based on a 26cm foot

Carbon Copy® High Performance » Toe Spring Stiffness Chart (K)

	Foot Size	Toe Spring Stiffness				
		1 (Very Low)	2 (Low)	3 (Regular)	4 (High)	5 (Very High)
Mens	23–24 cm	Less than 37kg (81lbs)	37–57 kg (81–125 lbs)	58–75 kg (126–165 lbs)	75–90 kg (166–200 lbs)	91–115kg (201–250 lbs)
	25–26 cm	Less than 42kg (91lbs)	42–59 kg (91–130 lbs)	60–77 kg (131–170 lbs)	78–93 kg (171–205 lbs)	94–115kg (206–250 lbs)
	27–28 cm	Less than 46kg (101lbs)	46–61 kg (101–135 lbs)	62–79 kg (136–175 lbs)	80–95 kg (176–210 lbs)	96–115kg (211–250 lbs)
	29–31 cm	Less than 51kg (111lbs)	51–66 kg (111–145 lbs)	67–82 kg (146–180 lbs)	83–98 kg (181–215 lbs)	99–115kg (216–250 lbs)
Womens	22–23 cm	Less than 33kg (71lbs)	33–52 kg (71–115 lbs)	53–73 kg (116–160 lbs)	74–88 kg (161–195 lbs)	89–115 kg (196–250 lbs)
	24–25 cm	Less than 37kg (81lbs)	37–57 kg (81–125 lbs)	58–75 kg (126–165 lbs)	75–90 kg (166–200 lbs)	91–115 kg (201–250 lbs)
	26–27 cm	Less than 42kg (91lbs)	42–59 kg (91–130 lbs)	60–77 kg (131–170 lbs)	78–93 kg (171–205 lbs)	94–115kg (206–250 lbs)

If the patient's weight is near the border between two categories:

a. choose the higher of the two categories if the patient is very active or has a **longer** residual limb

b. choose the lower of the two categories if the patient is moderately active or has a **very short** residual limb