

» DuraLite™ Foot

Prosthetic Foot

Supplier:

WillowWood™

Features & Benefits:

- Lightweight
- Unique composites provide amputees with increased deflection and energy return
- Field adjustable height to suit patient needs (13" maximum to 5.9" minimum)
- Adjustable heel spring options. Interchangeable heel wedges provided.
- Six toe resistance options
- Three proximal adapter options
- Urethane toe pad enhances medial and lateral motion
- Split heel conforms to the floor and provides stability on uneven terrain
- Stair Guard: Helps to reduce the risk of the shin adaptor catching on a step when descending stairs

Specifications:

- Types: *DuraLite foot with Pyramid adapter, Pyramid receiver or Torsion receiver*
- Sizes: 23–31 cm
- Weight limit: *Foot size 23–24 cm = 104kg (230lbs)*
Foot size 25–26 cm = 122kg (270lbs)
Foot size 27–31 cm = 136kg (300lbs)
- Weight: *356g (12.5oz) (Based on size 26cm foot, without accessories)*
- Toe spring stiffnesses: 1–6 (See chart)
- Heel heights: *10mm (3/8") or 19mm (3/4")*
- Clearances: *Maximum = 330mm (13")*
Pyramid adapter = 3/8" heel — 166mm (6 9/16")
3/4" heel — 156mm (6 1/8")
Pyramid receiver = 3/8" heel — 163mm (6 7/16")
3/4" heel — 152mm (6")
Torsion receiver = 3/8" heel — 208mm (8 3/16")
3/4" heel — 198mm (7 13/16")
- Heel insert densities: *Soft, Medium and Firm*
- Foot shell: *Lifelike in Buff, Tan* or Medium Brown* (* Special order only)*
- Activity level: *High (K3)*
- Includes foot shell, Spectra sock, choice of proximal adapter, 1 of 6 toe spring options and 3 wedges
- 46mm and 24mm height extension adapters available for use with torsion receiver

Indications:

- Transtibial and transfemoral amputees

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- 36 months on the foot, against manufacturers defects
- 9 months on foot shell, against manufacturers defects
- 3 weeks from date of purchase for installation



1 = DuraLite foot with Pyramid adapter

2 = Pyramid Receiver for Duralite foot



3 = DuraLite foot with torsion receiver

DuraLite™ Foot » WillowWood™

Ordering example: OW-DLF-24LM-3L-2

Item No.	Description	Side	Color	Heel Height	Weight ¹
OW-DLF-SSLB-KL-A	Prosthetic foot	Left	Buff	10mm (3/8")	356g (12.5oz)
OW-DLF-SSRB-KL-A	Prosthetic foot	Right	Buff	10mm (3/8")	356g (12.5oz)
OW-DLF-SSLM-KL-A	Prosthetic foot	Left	Medium brown ²	10mm (3/8")	356g (12.5oz)
OW-DLF-SSRM-KL-A	Prosthetic foot	Right	Medium brown ²	10mm (3/8")	356g (12.5oz)
OW-DLF-SSLT-KL-A	Prosthetic foot	Left	Tan ²	10mm (3/8")	356g (12.5oz)
OW-DLF-SSRT-KL-A	Prosthetic foot	Right	Tan ²	10mm (3/8")	356g (12.5oz)
OW-DLF-SSLB-KM-A	Prosthetic foot	Left	Buff	19mm (3/4")	356g (12.5oz)
OW-DLF-SSRB-KM-A	Prosthetic foot	Right	Buff	19mm (3/4")	356g (12.5oz)
OW-DLF-SSLM-KM-A	Prosthetic foot	Left	Medium brown ²	19mm (3/4")	356g (12.5oz)
OW-DLF-SSRM-KM-A	Prosthetic foot	Right	Medium brown ²	19mm (3/4")	356g (12.5oz)
OW-DLF-SSLT-KM-A	Prosthetic foot	Left	Tan ²	19mm (3/4")	356g (12.5oz)
OW-DLF-SSRT-KM-A	Prosthetic foot	Right	Tan ²	19mm (3/4")	356g (12.5oz)

SS = Size (23, 24, 25, 26, 27, 28, 29, 30, 31) **K** = Toe Spring (1–6) **A** = Adapter (1 [Pyramid Adapter], 2 [Pyramid Receiver], 3 [Torsion Receiver])

¹ Based on a 26cm foot

² Special order only

DuraLite™ » Composite Toe Spring Stiffness Chart (K)

Body Weight	Activity	Foot Size and Toe Spring Stiffness		
		23–24 cm	25–26 cm	27–31 cm
45–54 kg (100–120 lbs)	Medium	1	1	1
	High	2	2	2
55–68 kg (121–150 lbs)	Low	1	1	1
	Medium	2	2	2
	High	3	3	3
69–86 kg (151–190 lbs)	Low	2	2	2
	Medium	3	3	3
	High	4	4	4
87–104 kg (191–230 lbs)	Low	3	3	3
	Medium	4	4	4
	High	5	5	5
105–122 kg (231–270 lbs)	Low	—	4	4
	Medium	—	5	5
	High	—	—	6
123–135 kg (271–300 lbs)	Medium	—	—	5
	High	—	—	6

DuraLite™ » Height Extension Adapters

Item No.	Description	Material	Height	Weight Limit	Weight
OW-FND-136124	Height extension adapter	Titanium	24mm (15/16")	136kg (300lbs)	79g (2.8oz)
OW-FND-136146	Height extension adapter	Titanium	46mm (1 13/16")	136kg (300lbs)	96g (3.4oz)

Note: For use with the torsion receiver adapter only