

» SACH Foot

Prosthetic Foot

Supplier:

WillowWood™

Features & Benefits:

- Lightweight construction
- Attractive lifelike unisex cosmesis
- Water-resistant
- Lightweight composite keel
- Interchangeable heel inserts
- Flexible toe belting with three toe resistances

Specifications:

- Models: *SACH Foot (without pyramid)*
SACH Foot with Titanium Pyramid
- Sizes: 22–31 cm
- Weight limit: 115kg (250lbs)
- Weight: *SACH Foot* = 490g (17.3oz)
SACH Foot with Titanium Pyramid = 476g (17.1oz)
(Based on 26cm size foot)
- Toe spring stiffnesses: 2–4 (See chart)
- Heel heights: 10mm ($\frac{3}{8}$ ") or 19mm ($\frac{3}{4}$ ")
- Clearance: *SACH Foot* = 86mm ($3\frac{3}{8}$ ")
SACH Foot with Titanium Pyramid = 89mm ($3\frac{1}{2}$ ")
- Heel insert densities: *Soft, Medium and Firm*
- Foot shell: *Unisex, Lifelike in Buff, Tan or Medium Brown*
- Activity level: *Low (K1)*
- Includes all three heel inserts

Indications:

- Ideal for low activity patients

Contraindications:

- High activity level

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- One year and three weeks from purchase, against manufacturers defects



Without Pyramid



With Titanium Pyramid

SACH Foot » WillowWood™

Ordering example: OW-SFT-24LM-3L

Item No.	Description	Side	Color	Heel Height	Weight ¹
SACH Foot					
OW-SFT- SSLB-KL	Prosthetic foot	Left	Buff	10mm (3/8")	476g (16.8oz)
OW-SFT- SSRB-KL	Prosthetic foot	Right	Buff	10mm (3/8")	476g (16.8oz)
OW-SFT- SSLM-KL	Prosthetic foot	Left	Medium Brown	10mm (3/8")	476g (16.8oz)
OW-SFT- SSRM-KL	Prosthetic foot	Right	Medium Brown	10mm (3/8")	476g (16.8oz)
OW-SFT- SSLT-KL	Prosthetic foot	Left	Tan	10mm (3/8")	476g (16.8oz)
OW-SFT- SSRT-KL	Prosthetic foot	Right	Tan	10mm (3/8")	476g (16.8oz)
OW-SFT- SSLB-KM	Prosthetic foot	Left	Buff	19mm (3/4")	476g (16.8oz)
OW-SFT- SSRB-KM	Prosthetic foot	Right	Buff	19mm (3/4")	476g (16.8oz)
OW-SFT- SSLM-KM	Prosthetic foot	Left	Medium Brown	19mm (3/4")	476g (16.8oz)
OW-SFT- SSRM-KM	Prosthetic foot	Right	Medium Brown	19mm (3/4")	476g (16.8oz)
OW-SFT- SSLT-KM	Prosthetic foot	Left	Tan	19mm (3/4")	476g (16.8oz)
OW-SFT- SSRT-KM	Prosthetic foot	Right	Tan	19mm (3/4")	476g (16.8oz)

SACH Foot with titanium pyramid

OW-SFP- SSLB-KL	Prosthetic foot with titanium pyramid	Left	Buff	10mm (3/8")	490g (17.3oz)
OW-SFP- SSRB-KL	Prosthetic foot with titanium pyramid	Right	Buff	10mm (3/8")	490g (17.3oz)
OW-SFP- SSLM-KL	Prosthetic foot with titanium pyramid	Left	Medium Brown	10mm (3/8")	490g (17.3oz)
OW-SFP- SSRM-KL	Prosthetic foot with titanium pyramid	Right	Medium Brown	10mm (3/8")	490g (17.3oz)
OW-SFP- SSLT-KL	Prosthetic foot with titanium pyramid	Left	Tan	10mm (3/8")	490g (17.3oz)
OW-SFP- SSRT-KL	Prosthetic foot with titanium pyramid	Right	Tan	10mm (3/8")	490g (17.3oz)
OW-SFP- SSLB-KM	Prosthetic foot with titanium pyramid	Left	Buff	19mm (3/4")	490g (17.3oz)
OW-SFP- SSRB-KM	Prosthetic foot with titanium pyramid	Right	Buff	19mm (3/4")	490g (17.3oz)
OW-SFP- SSLM-KM	Prosthetic foot with titanium pyramid	Left	Medium Brown	19mm (3/4")	490g (17.3oz)
OW-SFP- SSRM-KM	Prosthetic foot with titanium pyramid	Right	Medium Brown	19mm (3/4")	490g (17.3oz)
OW-SFP- SSLT-KM	Prosthetic foot with titanium pyramid	Left	Tan	19mm (3/4")	490g (17.3oz)
OW-SFP- SSRT-KM	Prosthetic foot with titanium pyramid	Right	Tan	19mm (3/4")	490g (17.3oz)

SS= Size (22, 23, 24, 25, 26, 27, 28, 29, 30, 31)

K= Toe Spring (2–4)

¹ Based on a 26cm foot

SACH Foot » Toe Resistance Chart (K)

	2 Low	3 Regular	4 High
kg (lbs)	59kg & Below (130lbs & Below)	60–81 kg (131–180 lbs)	82–115 kg (181–250 lbs)