

EasyLiner

Liner

Supplier:

ALPS South Corporation

Features & Benefits:

- Reinforced seams improve durability and appearance
- Re-engineered distal end maximizes stability and strength without compromising comfort
- High performance fabric improves durability
- Longer lengths to accommodate more users

Specifications:

- Liner types available: *Cushion* or *Locking*
- Styles available: *Tapered* and *Uniform*
- Gel thickness available:
Uniform = 3 and 6 mm
Tapered = from 6 mm at distal end to 3 mm proximal end
- Sizes available: From 16 cm to 44 cm
- Colors available: Beige
- Measurements taken 6 cm from distal end

Indications:

- Transtibial amputations

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- 12 months on a single and 24 months on a pair, against manufacturers defects



768

Ordering Guidelines			
<i>Sample Order: ALPS-ELDT-16-3 ▼</i>			
ALPS	ELDT	16	3
	Model	Size	Gel Thickness
	ELDT = EasyLiner Locking Liner	16 = 16–19 cm	3 = 3mm (Uniform)
		20 = 20–23 cm	
		24 = 24–25 cm	6 = 6mm (Uniform)
		26 = 26–27 cm	
	ELFR = EasyLiner Cushion Liner	28 = 28–31 cm	6/3 = 6 to 3 mm (Tapered)
		32 = 32–35 cm	
		38 = 36–44 cm	
		44 = 43–60 cm	

EasyLiner » ALPS South Corporation

Style Sizing Chart				
Style	Gel Thickness	Size	*Distal End Circumference	Length
Tapered	6 to 3 mm	16	16–19 cm	33.0–35.6 cm
		20	20–23 cm	43.2–45.7 cm
		24	24–25 cm	41.9–45.7 cm
		26	26–27 cm	41.9–45.7 cm
		28	28–31 cm	43.2–45.7 cm
		32	32–35 cm	43.2–45.7 cm
		38	36–44 cm	43.2–45.7 cm
		44	43–60 cm	43.2–45.7 cm
Uniform	3mm, 6mm	16	16–19 cm	33.0–35.6 cm
		20	20–23 cm	43.2–45.7 cm
		24	24–25 cm	43.2–45.7 cm
		26	26–27 cm	43.2–45.7 cm
		28	28–31 cm	43.2–45.7 cm
		32	32–35 cm	38.1–45.7 cm
		38	36–44 cm	43.2–45.7 cm
		44	43–60 cm	43.2–45.7 cm

**Measure distal end circumference approximately 6 cm from distal stump end*