

## » Extreme Gel Liner

### Liner

#### Supplier:

ALPS South Corporation

#### Features & Benefits:

- Same gel properties as EasyGel
- Limited vertical stretch (80% less than other ALPS gel liners) reduces movement of redundant tissue
- High circumferential stretch to assist with form fitting
- Wider distal end for ease of donning
- Distal ends are proportionate to the size of the liner

#### Specifications:

- Types: *Cushion and Locking*
- Style: *Uniform*
- Material: *Proprietary gel*
- Thicknesses: *3mm and 6mm*
- Sizes: *16, 20, 24, 26, 28, 32, 38 and 44*
- Measurements taken 6cm from distal end

#### Indications:


- Transfemoral amputee
- Active transtibial amputee

#### Call Points:

- Prosthetist
- Prosthetic technician

#### Warranty:

- 12 months for single liner or 24 months for two identical liners, against manufacturers defects



**Extreme Gel Liner » ALPS South Corporation**

Item No.	Description	Style	Thickness	Size	Distal Circ*	Length
ALPS-AKDT-16-3	Locking gel liner	Uniform	3mm	16	16–19 cm	33–35.6 cm (13–14")
ALPS-AKDT-20-3	Locking gel liner	Uniform	3mm	20	20–23 cm	43.2–45.7 cm (17–18")
ALPS-AKDT-24-3	Locking gel liner	Uniform	3mm	24	24–25 cm	43.2–45.7 cm (17–18")
ALPS-AKDT-26-3	Locking gel liner	Uniform	3mm	26	26–27 cm	43.2–45.7 cm (17–18")
ALPS-AKDT-28-3	Locking gel liner	Uniform	3mm	28	28–31 cm	43.2–45.7 cm (17–18")
ALPS-AKDT-32-3	Locking gel liner	Uniform	3mm	32	32–37 cm	43.2–45.7 cm (17–18")
ALPS-AKDT-38-3	Locking gel liner	Uniform	3mm	38	38–43 cm	43.2–45.7 cm (17–18")
ALPS-AKDT-44-3	Locking gel liner	Uniform	3mm	44	44–53 cm	43.2–45.7 cm (17–18")
ALPS-AKDT-16-6	Locking gel liner	Uniform	6mm	16	16–19 cm	33–35.6 cm (13–14")
ALPS-AKDT-20-6	Locking gel liner	Uniform	6mm	20	20–23 cm	43.2–45.7 cm (17–18")
ALPS-AKDT-24-6	Locking gel liner	Uniform	6mm	24	24–25 cm	43.2–45.7 cm (17–18")
ALPS-AKDT-26-6	Locking gel liner	Uniform	6mm	26	26–27 cm	43.2–45.7 cm (17–18")
ALPS-AKDT-28-6	Locking gel liner	Uniform	6mm	28	28–31 cm	43.2–45.7 cm (17–18")
ALPS-AKDT-32-6	Locking gel liner	Uniform	6mm	32	32–37 cm	43.2–45.7 cm (17–18")
ALPS-AKDT-38-6	Locking gel liner	Uniform	6mm	38	38–43 cm	43.2–45.7 cm (17–18")
ALPS-AKDT-44-6	Locking gel liner	Uniform	6mm	44	44–53 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-16-3	Cushion gel liner	Uniform	3mm	16	16–19 cm	33–35.6 cm (13–14")
ALPS-AKFR-20-3	Cushion gel liner	Uniform	3mm	20	20–23 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-24-3	Cushion gel liner	Uniform	3mm	24	24–25 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-26-3	Cushion gel liner	Uniform	3mm	26	26–27 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-28-3	Cushion gel liner	Uniform	3mm	28	28–31 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-32-3	Cushion gel liner	Uniform	3mm	32	32–37 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-38-3	Cushion gel liner	Uniform	3mm	38	38–43 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-44-3	Cushion gel liner	Uniform	3mm	44	44–53 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-16-6	Cushion gel liner	Uniform	6mm	16	16–19 cm	33–35.6 cm (13–14")
ALPS-AKFR-20-6	Cushion gel liner	Uniform	6mm	20	20–23 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-24-6	Cushion gel liner	Uniform	6mm	24	24–25 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-26-6	Cushion gel liner	Uniform	6mm	26	26–27 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-28-6	Cushion gel liner	Uniform	6mm	28	28–31 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-32-6	Cushion gel liner	Uniform	6mm	32	32–37 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-38-6	Cushion gel liner	Uniform	6mm	38	38–43 cm	43.2–45.7 cm (17–18")
ALPS-AKFR-44-6	Cushion gel liner	Uniform	6mm	44	44–53 cm	43.2–45.7 cm (17–18")

\* Measure distal end circumference approximately 6cm from distal stump end