

ALPS Insoles

Foot Orthotics

Supplier:

ALPS South Corporation

Features & Benefits:

- Thermoformable
- Correction of foot pronation and supination
- Pressure distribution
- Treats transverse and longitudinal platypodia by modifying lower foam

Specifications:

- Available with or without metatarsal support
- Tri-laminate with rigid material located between two layers of foam
- Perforated top cover
- Durometer for top layer: Shore A 30
- Durometer for bottom layer: Shore A 50
- Heating temperature: 115°C (240°F)
- Heating time: 3 minutes

Indications:

- Treats transverse and longitudinal platypodia

Call Points:

- Pedorthist
- Orthotist
- Physiotherapist
- Occupational therapist

Warranty:

- Against manufacturers defects



Note: Thermoforming instructions for ALPS Insoles, are available

Item No.	Description
ALPS-ICNR-xx-2/3-HD	High Density Insole with 2/3 length foam and no metatarsal support
ALPS-ICMR-xx-2/3-HD	High Density Insole with 2/3 length foam and metatarsal support
ALPS-ICNR-xx-HD	Pediatric High Density Insole with full length foam and no metatarsal support (<i>only available in Size Code 23-40</i>)

xx indicates size (Size Code), see sizing chart on next page

ALPS Insoles » ALPS South Corporation

SIZING CHART			
Euro Size	Size		
(Size Code)	Men's	Women's	Youth
23	-	-	7
24	-	-	8
25	-	-	9.5
26	-	-	10.5
27	-	-	11.5
28	-	-	12.5
29	-	-	13
30	-	-	1
31	-	-	2
32	-	-	2.5
33	-	-	3
34	-	-	-
35	3.5	5	-
36	4.5	6	-
37	5	7	-
38	6	7.5	-
39	7	8	-
40	7.5	9	-
41	8.5	9.5	-
42	9	10	-
43	9.5	11	-
44	10	12	-
45	11	13	-
46	12	14	-

