

» Dynamic Walk / Dynamic Walk Single Side

Dynamic AFO

Supplier:

Centri AB (Fillauer)

Features & Benefits:

- Provides dorsi-assist with subtalar inversion and eversion
- Open heel construction allows natural ankle movement
- Excellent medial/lateral stability
- Flexible calf band

Specifications:

- Models: *Dynamic Walk*
Dynamic Walk Single Side
- Styles: *Standard= Regular and Wide ankle*
Single side= Regular
- Sizes: *Small, Medium and Large*
- Materials: *Comfil® Thermoformable Composite with PEEK rods*
(Note: PEEK rods on either the lateral or medial side are available for the Dynamic Walk Single Side)
- Thermoplastic carbon fiber can be reshaped by spot heating with heatgun (do not exceed 120°C)
- When selecting the appropriate size Dynamic Walk or Dynamic Walk Wide Ankle, the appropriate width measurement should be taken in the area between the PEEK insertion
- Do not heat or reshape PEEK rods or the surrounding area where they attach
- Available with custom sized calf bands (longer lead time required)

Indications:

- Peroneal paralysis
- Mild to moderate foot drop

Contraindications:

- Patients with contractures
- Ankle instability
- Patients with muscle spasticity

Call Points:

- Orthotist
- Physiotherapist
- Occupational therapist
- General practitioner
- Neurologist
- Neurosurgeon

Warranty:

- One year, against manufacturer defects



Dynamic Walk



Dynamic Walk Single Side

Dynamic Walk / Dynamic Walk Single Side » Centri AB (Fillauer)

Item No.	Description	Size	Side
Dynamic Walk — Standard			
DF-70011121	Dynamic dorsiflexion assist AFO without padding	Small	Left
DF-70011122	Dynamic dorsiflexion assist AFO without padding	Small	Right
DF-70011141	Dynamic dorsiflexion assist AFO without padding	Medium	Left
DF-70011142	Dynamic dorsiflexion assist AFO without padding	Medium	Right
DF-70011161	Dynamic dorsiflexion assist AFO without padding	Large	Left
DF-70011162	Dynamic dorsiflexion assist AFO without padding	Large	Right
DF-70011321	Wide ankle dynamic dorsiflexion assist AFO without padding	Small	Left
DF-70011322	Wide ankle dynamic dorsiflexion assist AFO without padding	Small	Right
DF-70011341	Wide ankle dynamic dorsiflexion assist AFO without padding	Medium	Left
DF-70011342	Wide ankle dynamic dorsiflexion assist AFO without padding	Medium	Right
DF-70011361	Wide ankle dynamic dorsiflexion assist AFO without padding	Large	Left
DF-70011362	Wide ankle dynamic dorsiflexion assist AFO without padding	Large	Right

Dynamic Walk — Single Side

DF-700113121	Medial side dynamic dorsiflexion assist AFO without padding	Small	Left
DF-700113122	Medial side dynamic dorsiflexion assist AFO without padding	Small	Right
DF-700113141	Medial side dynamic dorsiflexion assist AFO without padding	Medium	Left
DF-700113142	Medial side dynamic dorsiflexion assist AFO without padding	Medium	Right
DF-700113161	Medial side dynamic dorsiflexion assist AFO without padding	Large	Left
DF-700113162	Medial side dynamic dorsiflexion assist AFO without padding	Large	Right
DF-700114121	Lateral side dynamic dorsiflexion assist AFO without padding	Small	Left
DF-700114122	Lateral side dynamic dorsiflexion assist AFO without padding	Small	Right
DF-700114141	Lateral side dynamic dorsiflexion assist AFO without padding	Medium	Left
DF-700114142	Lateral side dynamic dorsiflexion assist AFO without padding	Medium	Right
DF-700114161	Lateral side dynamic dorsiflexion assist AFO without padding	Large	Left
DF-700114162	Lateral side dynamic dorsiflexion assist AFO without padding	Large	Right

Dynamic Walk — Custom Configuration

DF-700100000	Dynamic dorsiflexion assist AFO with custom calf band, without padding*	—	—
--------------	---	---	---

* When ordering custom configuration, reference a standard model item number with the required size for calf band








Dynamic Walk / Dynamic Walk Single Side » Padding

Item No.	Description	Size
DF-700100220	Padding for Dynamic Walk	Small
DF-700100240	Padding for Dynamic Walk	Medium
DF-700100260	Padding for Dynamic Walk	Large








Dynamic Walk / Dynamic Walk Single Side » Standard Sizing Information

Please note: The sizing information listed below are key in the successful fitting of the Dynamic Walk.








Dynamic Walk

							
Size	Foot Plate Length	Height	Width ¹	Calf Diameter	Size (Women's) ²	Size (Men's) ²	Weight
Small	230mm (9")	320mm (12 ¹⁹ / ₃₂ "	67.5mm (2 ²³ / ₃₂ "	82mm (3 ¹ / ₄ "	5½–8	4–6½	105g (3.7oz)
Medium	250mm (9 ²⁷ / ₃₂ "	340mm (13 ¹³ / ₃₂ "	72.5mm (2 ²⁷ / ₃₂ "	100mm (3 ³¹ / ₃₂ "	8½–10½	7–9	122g (4.3oz)
Large	270mm (9 ⁵ / ₈ "	370mm (14 ⁹ / ₁₆ "	80.5mm (3 ⁵ / ₃₂ "	120mm (4 ²³ / ₃₂ "	12 – 15½	10½–14	140g (4.9oz)

Dynamic Walk Wide Ankle

							
Size	Foot Plate Length	Height	Width ¹	Calf Diameter	Size (Women's) ²	Size (Men's) ²	Weight
Small	230mm (9")	320mm (12 ¹⁹ / ₃₂ "	72.5mm (2 ²⁷ / ₃₂ "	82mm (3 ¹ / ₄ "	5½–8	4–6½	106g (3.7oz)
Medium	250mm (9 ²⁷ / ₃₂ "	340mm (13 ¹³ / ₃₂ "	80.5mm (3 ⁵ / ₃₂ "	100mm (3 ³¹ / ₃₂ "	8½–10½	7–9	123g (4.3oz)
Large	270mm (9 ⁵ / ₈ "	370mm (14 ⁹ / ₁₆ "	88.5mm (3 ¹ / ₂ "	120mm (4 ²³ / ₃₂ "	12 – 15½	10½–14	142g (5oz)

Dynamic Walk Single Side

							
Size	Foot Plate Length	Height	Width ¹	Calf Diameter	Size (Women's) ²	Size (Men's) ²	Weight
Small	230mm (9")	340mm (13 ¹³ / ₃₂ "	—	82mm (3 ¹ / ₄ "	5½–8	4–6½	103g (3.6oz)
Medium	250mm (9 ²⁷ / ₃₂ "	370mm (14 ⁹ / ₁₆ "	—	100mm (3 ³¹ / ₃₂ "	8½–10½	7–9	120g (4.2oz)
Large	270mm (9 ⁵ / ₈ "	385mm (15 ⁵ / ₁₆ "	—	120mm (4 ²³ / ₃₂ "	12 – 15½	10½–14	140g (4.9oz)

¹ = Width of area between PEEK insertion

² = Approximate sizes only. For proper fitting, consult with patient or measure their shoe size prior to ordering.