

» Element Foot System™

Prosthetic Foot

Supplier:

Fillauer LLC

Features & Benefits:

- Unique 3rd Element technology allows for a more natural gait cycle
- Smooth rollover and superior stability
- Interchangeable heel elastomers for fine tuning heel stiffness
- Modular design with superior flexibility
- Durable and lightweight
- Innovative foot shell — flexible with realistic features

Specifications:

- Sizes: 22–30 cm
- Weight limit: 125kg (275lbs)
- Weight: 340g (12oz)
(Based on size 26cm foot, without foot shell)
- Spring stiffnesses: 11, 22, 32, 43, 53, 54, 64, 65 and 75 (see chart)
- Heel height: 10mm ($\frac{3}{8}$ ")
- Clearance: 95mm ($3\frac{3}{4}$ ")
- Foot shell: *Caucasian* (custom colors are available for order)
- Activity level: *High (K3)*
- Complies with ISO-22675 test standard
- Includes heel bands, Spectra sock and instruction manual

Indications:

- Transfemoral and transtibial amputees
- Active K3 amputees

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- 36 months on Element foot, against manufacturers defects
- Six months on foot shell, against manufacturers defects



Element Foot System™ » Fillauer LLC

Ordering example: DF-100-10-24-A22

Item No.	Description	Size	Clearance	Heel Height	Weight ¹
DF-100-10-22- KK	Low profile prosthetic foot	22cm	95mm (3¾")	10mm (¾")	340g (12oz)
DF-100-10-23- KK	Low profile prosthetic foot	23cm	95mm (3¾")	10mm (¾")	340g (12oz)
DF-100-10-24- KK	Low profile prosthetic foot	24cm	95mm (3¾")	10mm (¾")	340g (12oz)
DF-100-10-25- KK	Low profile prosthetic foot	25cm	95mm (3¾")	10mm (¾")	340g (12oz)
DF-100-10-26- KK	Low profile prosthetic foot	26cm	95mm (3¾")	10mm (¾")	340g (12oz)
DF-100-10-27- KK	Low profile prosthetic foot	27cm	95mm (3¾")	10mm (¾")	340g (12oz)
DF-100-10-28- KK	Low profile prosthetic foot	28cm	95mm (3¾")	10mm (¾")	340g (12oz)
DF-100-10-29- KK	Low profile prosthetic foot	29cm	95mm (3¾")	10mm (¾")	340g (12oz)
DF-100-10-30- KK	Low profile prosthetic foot	30cm	95mm (3¾")	10mm (¾")	340g (12oz)

KK = Spring. See the below chart for spring selection.

¹ Based on a 26cm foot, without foot shell

Element » Heel Spring Stiffness (KK)

	45–54 kg (100–119 lbs)	55–63 kg (120–139 lbs)	64–72 kg (140–159 lbs)	73–81 kg (160–179 lbs)	82–90 kg (180–199 lbs)	91–101 kg (200–224 lbs)	102–113 kg (225–249 lbs)	114–125 kg (250–275 lbs)
22cm	A11	A11	A11	A22	A22	A22	—	—
23cm	A11	A11	A11	A22	A22	A22	—	—
24cm	A11	A11	A11	A22	A22	A22	—	—
25cm	A32	A32	A32	A43	B43	B53	C54	C64
26cm	A32	A32	A32	A43	B43	B53	C54	C64
27cm	*	A43	A43	A43	B54	B54	C65	C65
28cm	*	A43	A43	A43	B54	B54	C65	C65
29cm	*	*	C64	C64	C64	D75	D75	D75
30cm	*	*	C64	C64	C64	D75	D75	D75

* Please contact your Inside Sales Representative for recommendation

Centri® MCV Foot Shell

DF-XXXX13033

	22cm	23cm	24cm	25cm	26cm	27cm	28cm	29cm	30cm
Left	4522	4523	4524	4525	4526	4527	4528	4529	4530
Right	4622	4623	4624	4625	4626	4627	4628	4629	4630

Note: Please allow for a 3 week lead time with custom colored foot shell orders