

» Element^{DS} Foot System™

Prosthetic Foot

Supplier:

Fillauer LLC

Features & Benefits:

- Unique 3rd Element technology allows for a more natural gait cycle
- Smooth rollover and superior stability
- Interchangeable heel elastomers for fine tuning heel stiffness
- Modular design with superior flexibility
- Durable and lightweight
- Vertical shock and 60° of axial rotation
- Innovative foot shell — flexible with realistic features

Specifications:

- Sizes: 22–30 cm
- Weight limit: 125kg (275lbs)
- Weight: 632g (22.3oz)
(Based on size 26cm foot, without foot shell)
- Spring stiffnesses: 11, 22, 32, 43, 53, 54, 64, 65 and 75
(see chart for more information on spring/shock combinations)
- Heel height: 10mm (³/₈"
- Clearance: 165mm (6½")
- Foot shell: *Caucasian* (custom colors are available for order)
- Activity level: *High (K3)*
- Complies with ISO-22675 test standard
- Includes Durashock pylon, heel bands and elastomers, Spectra sock and instruction manual

Indications:

- Transfemoral and transtibial amputees
- Active K3 amputees

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- 36 months on foot, against manufacturers defects
- Six months on foot shell, against manufacturers defects



Element^{DS} Foot System™ » Fillauer LLC

Ordering example: DF-105-10-24-A22B

Item No.	Description	Size	Clearance	Heel Height	Weight ¹
DF-105-10-22- KK	Low profile prosthetic foot with pylon	22cm	165mm (6½")	10mm (¾")	632g (22.3oz)
DF-105-10-23- KK	Low profile prosthetic foot with pylon	23cm	165mm (6½")	10mm (¾")	632g (22.3oz)
DF-105-10-24- KK	Low profile prosthetic foot with pylon	24cm	165mm (6½")	10mm (¾")	632g (22.3oz)
DF-105-10-25- KK	Low profile prosthetic foot with pylon	25cm	165mm (6½")	10mm (¾")	632g (22.3oz)
DF-105-10-26- KK	Low profile prosthetic foot with pylon	26cm	165mm (6½")	10mm (¾")	632g (22.3oz)
DF-105-10-27- KK	Low profile prosthetic foot with pylon	27cm	165mm (6½")	10mm (¾")	632g (22.3oz)
DF-105-10-28- KK	Low profile prosthetic foot with pylon	28cm	165mm (6½")	10mm (¾")	632g (22.3oz)
DF-105-10-29- KK	Low profile prosthetic foot with pylon	29cm	165mm (6½")	10mm (¾")	632g (22.3oz)
DF-105-10-30- KK	Low profile prosthetic foot with pylon	30cm	165mm (6½")	10mm (¾")	632g (22.3oz)

KK = Spring / shock combination. See the below chart for spring / shock selections.

¹ Based on a 26cm foot, without foot shell

Element^{DS} » Heel Spring Stiffness / Shock Combination (KK)

	45–54 kg (100–119 lbs)	55–63 kg (120–139 lbs)	64–72 kg (140–159 lbs)	73–81 kg (160–179 lbs)	82–90 kg (180–199 lbs)	91–101 kg (200–224 lbs)	102–113 kg (225–249 lbs)	114–125 kg (250–275 lbs)
22cm	A11R	A11R	A11R	A22R	A22B	A22B	—	—
23cm	A11R	A11R	A11R	A22B	A22B	A22K	—	—
24cm	A11R	A11R	A11R	A22B	A22B	A22K	—	—
25cm	A32R	A32B	A32B	A43B	B43B	B53K	C54K	C64K
26cm	A32R	A32B	A32B	A43B	B43K	B53K	C54K	C64K
27cm	*	A43B	A43B	A43K	B54K	B54K	C65K	C65K
28cm	*	A43B	A43K	A43K	B54K	B54K	C65K	C65K
29cm	*	*	C64K	C64K	C64K	D75K	D75K	D75K
30cm	*	*	C64K	C64K	C64K	D75K	D75K	D75K

* Please contact your Inside Sales Representative for recommendation

Centri[®] MCV Foot Shell

DF-XXXX13033

	22cm	23cm	24cm	25cm	26cm	27cm	28cm	29cm	30cm
Left	4522	4523	4524	4525	4526	4527	4528	4529	4530
Right	4622	4623	4624	4625	4626	4627	4628	4629	4630

Note: Please allow for a 3 week lead time with custom colored foot shell orders