

» Wave Sport Foot System™

Prosthetic Foot

Supplier:

Fillauer LLC

Features & Benefits:

- Unique Wave-Spring technology provides shock absorption and stability
- "C-shaped" pylon provides greater dynamic response
- Full heel plate allows patient to use foot in foot shell and/or wear a standard athletic shoe
- Lightweight and durable carbon fiber construction
- Custom built for each patient
- Modular adapters available for easy alignment
- Innovative foot shell — flexible with realistic features

Specifications:

- Types: *Standard and Low profile*
- Sizes: *22–30 cm*
- Weight limit: *125kg (275lbs)*
- Weight: *550g (19.4oz) (Based on size 27cm foot, without foot shell)*
- Heel height: *10mm (3/8")*
- Clearance: *Standard= 240mm (9 1/2")*
Low profile= 200mm (8")
- Cosmesis: *Caucasian (custom colors are available for order)*
- Activity level: *High to Very high (K3–K4)*
- Complies with ISO-22675 test standard
- Please allow for a seven day lead time for item
- Optional pyramid adapters, Spectra sock and foot shell available (sold separately)

Indications:

- Active K3–K4 transtibial amputees

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- 36 months on foot, against manufacturers defects
- Six months on foot shell, against manufacturers defects



Wave Sport Foot System™ » Fillauer LLC

Please note: When ordering this item, it is important to indicate patients activity level, weight and foot size

Ordering example: DF-127-10-SH-CSF-25

Item No.	Description	Size	Clearance	Heel Height	Weight ¹
DF-127-10-SH-CSF-22	High-impact prosthetic foot	22cm	240mm (9½")	10mm (¾")	550g (19.4oz)
DF-127-10-SH-CSF-23	High-impact prosthetic foot	23cm	240mm (9½")	10mm (¾")	550g (19.4oz)
DF-127-10-SH-CSF-24	High-impact prosthetic foot	24cm	240mm (9½")	10mm (¾")	550g (19.4oz)
DF-127-10-SH-CSF-25	High-impact prosthetic foot	25cm	240mm (9½")	10mm (¾")	550g (19.4oz)
DF-127-10-SH-CSF-26	High-impact prosthetic foot	26cm	240mm (9½")	10mm (¾")	550g (19.4oz)
DF-127-10-SH-CSF-27	High-impact prosthetic foot	27cm	240mm (9½")	10mm (¾")	550g (19.4oz)
DF-127-10-SH-CSF-28	High-impact prosthetic foot	28cm	240mm (9½")	10mm (¾")	550g (19.4oz)
DF-127-10-SH-CSF-29	High-impact prosthetic foot	29cm	240mm (9½")	10mm (¾")	550g (19.4oz)
DF-127-10-SH-CSF-30	High-impact prosthetic foot	30cm	240mm (9½")	10mm (¾")	550g (19.4oz)
DF-127-10-LP-CSF-22	Low profile high-impact prosthetic foot	22cm	200mm (8")	10mm (¾")	550g (19.4oz)
DF-127-10-LP-CSF-23	Low profile high-impact prosthetic foot	23cm	200mm (8")	10mm (¾")	550g (19.4oz)
DF-127-10-LP-CSF-24	Low profile high-impact prosthetic foot	24cm	200mm (8")	10mm (¾")	550g (19.4oz)
DF-127-10-LP-CSF-25	Low profile high-impact prosthetic foot	25cm	200mm (8")	10mm (¾")	550g (19.4oz)
DF-127-10-LP-CSF-26	Low profile high-impact prosthetic foot	26cm	200mm (8")	10mm (¾")	550g (19.4oz)
DF-127-10-LP-CSF-27	Low profile high-impact prosthetic foot	27cm	200mm (8")	10mm (¾")	550g (19.4oz)
DF-127-10-LP-CSF-28	Low profile high-impact prosthetic foot	28cm	200mm (8")	10mm (¾")	550g (19.4oz)
DF-127-10-LP-CSF-29	Low profile high-impact prosthetic foot	29cm	200mm (8")	10mm (¾")	550g (19.4oz)
DF-127-10-LP-CSF-30	Low profile high-impact prosthetic foot	30cm	200mm (8")	10mm (¾")	550g (19.4oz)

¹ Based on a 27cm foot, without foot shell.

NOTE: Please allow for a seven day lead time for item.

Wave » Accessories (Sold Separately)

Item No.	Description	Material	Weight Limit	Weight
Pyramid Adapters				
HD-60254	4-hole pyramid with hole	Titanium	136kg (300lbs)	50g (1.8oz)
HD-60293	4-hole receiver	Titanium	136kg (300lbs)	82g (2.9oz)
Spectra Sock				
DF-192500	Spectra Sock			

Centri® MCV Foot Shell

DF-XXXX13033									
Side	22cm	23cm	24cm	25cm	26cm	27cm	28cm	29cm	30cm
Left	4522	4523	4524	4525	4526	4527	4528	4529	4530
Right	4622	4623	4624	4625	4626	4627	4628	4629	4630

Note: Please allow for a 3 week lead time with custom colored foot shell orders