

» Therafirm® Compression Socks

Support Socks

Supplier:

Knit-Rite, Inc.

Features & Benefits:

- Gradient compression
- Greatest amount of pressure is concentrated at the ankle and decreases gradually the length of the leg
- Wide range of compression for an array of various extremity problems
- Comfortable to wear
- Non-slip and non-restrictive knit-in band for greater comfort
- Available in an array of colors

Specifications:

- Types: *Therafirm®* and *Therafirm Light™*
- Sizes: *Small, Medium, Large and X-Large**
(*Available for Therafirm model only)
- Materials: *Therafirm®* = Two way stretch nylon and spandex
Therafirm Light™ = Lycra® and nylon blend
- Compression: *Therafirm®* = Mild (15–20 mmHg)
Moderate (20–30 mmHg)**
Firm (30–40 mmHg)**
(**Available for men's sizes only)
Therafirm Light™ = 10–15 mmHg
- Available in men's and women's sizes

Indications:

- Tired, aching and fatigued legs
- Improvement of blood circulation
- Various forms of ankle swelling
- Various degrees of varicosities
- Post-schlerotherapy applications
- Lymphedema

Call Points:

- Podiatrist
- General Practitioner
- Pedorthist
- Chiropracist

Warranty:

- Against manufacturers defects



Men's trouser socks



Women's trouser socks

Therafirm® Compression Socks » Knit-Rite, Inc.

| Item No. | Description | Compression | Style | Color | Size |
|-------------------|---------------------------------------|-------------|---------|--------------|---------|
| Therafirm® | | | | | |
| KR-68310 | Mild compression support trouser sock | 15–20 mmHg | Men's | Black | Small |
| KR-68311 | Mild compression support trouser sock | 15–20 mmHg | Men's | Black | Medium |
| KR-68312 | Mild compression support trouser sock | 15–20 mmHg | Men's | Black | Large |
| KR-68340 | Mild compression support trouser sock | 15–20 mmHg | Men's | Black | X-Large |
| KR-68364 | Mild compression support trouser sock | 15–20 mmHg | Men's | Charcoal | Small |
| KR-68365 | Mild compression support trouser sock | 15–20 mmHg | Men's | Charcoal | Medium |
| KR-68366 | Mild compression support trouser sock | 15–20 mmHg | Men's | Charcoal | Large |
| KR-68367 | Mild compression support trouser sock | 15–20 mmHg | Men's | Charcoal | X-Large |
| KR-68350 | Mild compression support trouser sock | 15–20 mmHg | Men's | Navy | Small |
| KR-68319 | Mild compression support trouser sock | 15–20 mmHg | Men's | Navy | Medium |
| KR-68320 | Mild compression support trouser sock | 15–20 mmHg | Men's | Navy | Large |
| KR-68341 | Mild compression support trouser sock | 15–20 mmHg | Men's | Navy | X-Large |
| KR-68353 | Mild compression support trouser sock | 15–20 mmHg | Men's | Brown | Small |
| KR-68313 | Mild compression support trouser sock | 15–20 mmHg | Men's | Brown | Medium |
| KR-68314 | Mild compression support trouser sock | 15–20 mmHg | Men's | Brown | Large |
| KR-68342 | Mild compression support trouser sock | 15–20 mmHg | Men's | Brown | X-Large |
| KR-68363 | Mild compression support trouser sock | 15–20 mmHg | Men's | Khaki | Small |
| KR-68326 | Mild compression support trouser sock | 15–20 mmHg | Men's | Khaki | Medium |
| KR-68327 | Mild compression support trouser sock | 15–20 mmHg | Men's | Khaki | Large |
| KR-68362 | Mild compression support trouser sock | 15–20 mmHg | Men's | Khaki | X-Large |
| KR-68360 | Mild compression support trouser sock | 15–20 mmHg | Men's | White | Small |
| KR-68324 | Mild compression support trouser sock | 15–20 mmHg | Men's | White | Medium |
| KR-68325 | Mild compression support trouser sock | 15–20 mmHg | Men's | White | Large |
| KR-68361 | Mild compression support trouser sock | 15–20 mmHg | Men's | White | X-Large |
| KR-68500 | Mild compression support trouser sock | 15–20 mmHg | Women's | Black | Small |
| KR-68501 | Mild compression support trouser sock | 15–20 mmHg | Women's | Black | Medium |
| KR-68502 | Mild compression support trouser sock | 15–20 mmHg | Women's | Black | Large |
| KR-68503 | Mild compression support trouser sock | 15–20 mmHg | Women's | Black | X-Large |
| KR-68505 | Mild compression support trouser sock | 15–20 mmHg | Women's | Brown | Small |
| KR-68506 | Mild compression support trouser sock | 15–20 mmHg | Women's | Brown | Medium |
| KR-68507 | Mild compression support trouser sock | 15–20 mmHg | Women's | Brown | Large |
| KR-68508 | Mild compression support trouser sock | 15–20 mmHg | Women's | Brown | X-Large |
| KR-68510 | Mild compression support trouser sock | 15–20 mmHg | Women's | Navy | Small |
| KR-68511 | Mild compression support trouser sock | 15–20 mmHg | Women's | Navy | Medium |
| KR-68512 | Mild compression support trouser sock | 15–20 mmHg | Women's | Navy | Large |
| KR-68513 | Mild compression support trouser sock | 15–20 mmHg | Women's | Navy | X-Large |
| KR-68525 | Mild compression support trouser sock | 15–20 mmHg | Women's | Khaki | Small |
| KR-68526 | Mild compression support trouser sock | 15–20 mmHg | Women's | Khaki | Medium |
| KR-68527 | Mild compression support trouser sock | 15–20 mmHg | Women's | Khaki | Large |
| KR-68528 | Mild compression support trouser sock | 15–20 mmHg | Women's | Khaki | X-Large |
| KR-68530 | Mild compression support trouser sock | 15–20 mmHg | Women's | Winter White | Small |
| KR-68531 | Mild compression support trouser sock | 15–20 mmHg | Women's | Winter White | Medium |
| KR-68532 | Mild compression support trouser sock | 15–20 mmHg | Women's | Winter White | Large |
| KR-68533 | Mild compression support trouser sock | 15–20 mmHg | Women's | Winter White | X-Large |
| KR-68515 | Mild compression support trouser sock | 15–20 mmHg | Women's | White | Small |
| KR-68516 | Mild compression support trouser sock | 15–20 mmHg | Women's | White | Medium |
| KR-68517 | Mild compression support trouser sock | 15–20 mmHg | Women's | White | Large |
| KR-68518 | Mild compression support trouser sock | 15–20 mmHg | Women's | White | X-Large |

Therafirm® Compression Socks » Knit-Rite, Inc.

| Item No. | Description | Compression | Style | Color | Size |
|----------|---|-------------|-------|----------|---------|
| KR-68315 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Black | Small |
| KR-68316 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Black | Medium |
| KR-68317 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Black | Large |
| KR-68318 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Black | X-Large |
| KR-68368 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Charcoal | Small |
| KR-68369 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Charcoal | Medium |
| KR-68370 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Charcoal | Large |
| KR-68371 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Charcoal | X-Large |
| KR-68346 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Navy | Small |
| KR-68303 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Navy | Medium |
| KR-68304 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Navy | Large |
| KR-68309 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Navy | X-Large |
| KR-68347 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Brown | Small |
| KR-68301 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Brown | Medium |
| KR-68302 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Brown | Large |
| KR-68328 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Brown | X-Large |
| KR-68349 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Khaki | Small |
| KR-68321 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Khaki | Medium |
| KR-68322 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Khaki | Large |
| KR-68323 | Moderate compression support trouser sock | 20–30 mmHg | Men's | Khaki | X-Large |
| KR-68331 | Moderate compression support trouser sock | 20–30 mmHg | Men's | White | Small |
| KR-68307 | Moderate compression support trouser sock | 20–30 mmHg | Men's | White | Medium |
| KR-68308 | Moderate compression support trouser sock | 20–30 mmHg | Men's | White | Large |
| KR-68330 | Moderate compression support trouser sock | 20–30 mmHg | Men's | White | X-Large |
| KR-69315 | Firm compression support trouser sock | 30–40 mmHg | Men's | Black | Small |
| KR-69316 | Firm compression support trouser sock | 30–40 mmHg | Men's | Black | Medium |
| KR-69317 | Firm compression support trouser sock | 30–40 mmHg | Men's | Black | Large |
| KR-69318 | Firm compression support trouser sock | 30–40 mmHg | Men's | Black | X-Large |
| KR-69346 | Firm compression support trouser sock | 30–40 mmHg | Men's | Navy | Small |
| KR-69303 | Firm compression support trouser sock | 30–40 mmHg | Men's | Navy | Medium |
| KR-69304 | Firm compression support trouser sock | 30–40 mmHg | Men's | Navy | Large |
| KR-69309 | Firm compression support trouser sock | 30–40 mmHg | Men's | Navy | X-Large |
| KR-69331 | Firm compression support trouser sock | 30–40 mmHg | Men's | White | Small |
| KR-69307 | Firm compression support trouser sock | 30–40 mmHg | Men's | White | Medium |
| KR-69308 | Firm compression support trouser sock | 30–40 mmHg | Men's | White | Large |
| KR-69330 | Firm compression support trouser sock | 30–40 mmHg | Men's | White | X-Large |

Therafirm® » Sizing Information

| Size | Foot Length | Calf Circumference | Ankle Circumference | Shoe Size |
|----------------|--------------------|----------------------|---------------------|-----------|
| Men's | | | | |
| Small | 28–36 cm (11"–14") | 27–38 cm (10½"–15") | 18–23 cm (7"–9") | 7 or less |
| Medium | 31–38 cm (12"–15") | 28–42 cm (11"–16½") | 23–28 cm (9"–11") | 7½–10 |
| Large | 33–41 cm (13"–16") | 31–43 cm (12"–17") | 25–36 cm (10"–14") | 10½–12 |
| X-Large | 36–46 cm (14"–18") | 33–48 cm (13"–19") | 28–38 cm (11"–15") | 12½ + |
| Women's | | | | |
| Small | 31–46 cm (12"–18") | 28–34 cm (11"–13½") | 16–18 cm (6¼"–7¼") | 4–5 |
| Medium | 33–48 cm (13"–19") | 33–39 cm (13"–15½") | 19–22 cm (7½"–9") | 5½–7½ |
| Large | 36–51 cm (14"–20") | 38–44 cm (15"–17½") | 24–27 cm (9¼"–10¾") | 8–10½ |
| X-Large | 38–56 cm (15"–22") | 45–47 cm (17¾"–22½") | 28–29 cm (11"–11½") | 11–12 |

Therafirm® Compression Socks » Knit-Rite, Inc.

| Item No. | Description | Compression | Style | Color | Size |
|-------------------------|--|-------------|---------|--------------|--------|
| Therafirm Light™ | | | | | |
| KR-90415 | Light compression support trouser sock | 10–15 mmHg | Men's | Black | Small |
| KR-90425 | Light compression support trouser sock | 10–15 mmHg | Men's | Black | Medium |
| KR-90435 | Light compression support trouser sock | 10–15 mmHg | Men's | Black | Large |
| KR-90401 | Light compression support trouser sock | 10–15 mmHg | Men's | Brown | Small |
| KR-90421 | Light compression support trouser sock | 10–15 mmHg | Men's | Brown | Medium |
| KR-90431 | Light compression support trouser sock | 10–15 mmHg | Men's | Brown | Large |
| KR-90417 | Light compression support trouser sock | 10–15 mmHg | Men's | Navy | Small |
| KR-90427 | Light compression support trouser sock | 10–15 mmHg | Men's | Navy | Medium |
| KR-90437 | Light compression support trouser sock | 10–15 mmHg | Men's | Navy | Large |
| KR-90215 | Ribbed light compression support trouser sock | 10–15 mmHg | Women's | Black | Small |
| KR-90225 | Ribbed light compression support trouser sock | 10–15 mmHg | Women's | Black | Medium |
| KR-90235 | Ribbed light compression support trouser sock | 10–15 mmHg | Women's | Black | Large |
| KR-90216 | Ribbed light compression support trouser sock | 10–15 mmHg | Women's | Cocoa | Small |
| KR-90226 | Ribbed light compression support trouser sock | 10–15 mmHg | Women's | Cocoa | Medium |
| KR-90236 | Ribbed light compression support trouser sock | 10–15 mmHg | Women's | Cocoa | Large |
| KR-90218 | Ribbed light compression support trouser sock | 10–15 mmHg | Women's | Coal | Small |
| KR-90228 | Ribbed light compression support trouser sock | 10–15 mmHg | Women's | Coal | Medium |
| KR-90238 | Ribbed light compression support trouser sock | 10–15 mmHg | Women's | Coal | Large |
| KR-90213 | Ribbed light compression support trouser sock | 10–15 mmHg | Women's | White | Small |
| KR-90223 | Ribbed light compression support trouser sock | 10–15 mmHg | Women's | White | Medium |
| KR-90233 | Ribbed light compression support trouser sock | 10–15 mmHg | Women's | White | Large |
| KR-95315 | Diamond light compression support trouser sock | 10–15 mmHg | Women's | Black | Small |
| KR-95325 | Diamond light compression support trouser sock | 10–15 mmHg | Women's | Black | Medium |
| KR-95335 | Diamond light compression support trouser sock | 10–15 mmHg | Women's | Black | Large |
| KR-95316 | Diamond light compression support trouser sock | 10–15 mmHg | Women's | Cocoa | Small |
| KR-95326 | Diamond light compression support trouser sock | 10–15 mmHg | Women's | Cocoa | Medium |
| KR-95336 | Diamond light compression support trouser sock | 10–15 mmHg | Women's | Cocoa | Large |
| KR-95317 | Diamond light compression support trouser sock | 10–15 mmHg | Women's | Navy | Small |
| KR-95327 | Diamond light compression support trouser sock | 10–15 mmHg | Women's | Navy | Medium |
| KR-95337 | Diamond light compression support trouser sock | 10–15 mmHg | Women's | Navy | Large |
| KR-95314 | Diamond light compression support trouser sock | 10–15 mmHg | Women's | Winter White | Small |
| KR-95324 | Diamond light compression support trouser sock | 10–15 mmHg | Women's | Winter White | Medium |
| KR-95334 | Diamond light compression support trouser sock | 10–15 mmHg | Women's | Winter White | Large |

Therafirm Light™ » Sizing Information

| Size | Shoe Size (Women's) | Shoe Size (Men's) |
|--------|---------------------|-------------------|
| Small | 5–8 | 7½–10 |
| Medium | 8½–11 | 10½–12 |
| Large | 7–10 | 12½ + |

*Fits up to 8½" calf circumference