

» Cellacast® Active

Semi-rigid Polyester Casting Impression

Supplier:

Lohmann & Rauscher, Inc.

Features & Benefits:

- Allows controlled movement while focused area is immobilized
- Applied as easily as standard Plaster of Paris bandages
- Thin material is highly conformable to body contours
- X-ray transparent
- Easily removed or worked on with cast scissors

Specifications:

- Materials: *Polyester fabric impregnated with polyurethane resin*
- Use of padding on patient is required prior to use
- Setting: *Using clean water (20–24°C), dip the bandage diagonally and lightly squeeze 2–3 times*
- Molding time: *2½–3 minutes*
- Weight bearing set time: *4–5 minutes*
- Full set time: *20–30 minutes*
- Storage: *Even temperature, no higher than 24°C*

Indications:


- Immobilization in primary or secondary treatment after fractures, operations, tendon injuries and overload conditions
- Partial immobilization
- Pediatric casts

Contraindications:

- Using water warmer than 24°C
- Avoid direct contact during bandaging—please use proper protection, in a well ventilated area

Warranty:

- Against manufacturers defects



Cellacast® Active » Lohmann & Rauscher, Inc.

| Item No. | Description | Color | Width | Length | Unit |
|-------------|--------------------------------------|--------|-------------|--------------|----------------|
| LR-23034-BX | Semi-rigid synthetic casting bandage | Creme | 2.5cm (1") | 1.8m (2 yds) | 10 Rolls / Box |
| LR-23035-BX | Semi-rigid synthetic casting bandage | Creme | 5cm (2") | 3.6m (4 yds) | 10 Rolls / Box |
| LR-23036-BX | Semi-rigid synthetic casting bandage | Creme | 7.5cm (3") | 3.6m (4 yds) | 10 Rolls / Box |
| LR-23037-BX | Semi-rigid synthetic casting bandage | Creme | 10cm (4") | 3.6m (4 yds) | 10 Rolls / Box |
| LR-23038-BX | Semi-rigid synthetic casting bandage | Creme | 12.5cm (5") | 3.6m (4 yds) | 10 Rolls / Box |
| LR-23040-BX | Semi-rigid synthetic casting bandage | Blue | 5cm (2") | 3.6m (4 yds) | 10 Rolls / Box |
| LR-23041-BX | Semi-rigid synthetic casting bandage | Blue | 7.5cm (3") | 3.6m (4 yds) | 10 Rolls / Box |
| LR-23042-BX | Semi-rigid synthetic casting bandage | Blue | 10cm (4") | 3.6m (4 yds) | 10 Rolls / Box |
| LR-23043-BX | Semi-rigid synthetic casting bandage | Blue | 12.5cm (5") | 3.6m (4 yds) | 10 Rolls / Box |
| LR-23045-BX | Semi-rigid synthetic casting bandage | Orange | 5cm (2") | 3.6m (4 yds) | 10 Rolls / Box |
| LR-23046-BX | Semi-rigid synthetic casting bandage | Orange | 7.5cm (3") | 3.6m (4 yds) | 10 Rolls / Box |
| LR-23047-BX | Semi-rigid synthetic casting bandage | Orange | 10cm (4") | 3.6m (4 yds) | 10 Rolls / Box |
| LR-23048-BX | Semi-rigid synthetic casting bandage | Orange | 12.5cm (5") | 3.6m (4 yds) | 10 Rolls / Box |