

## » Cellona® Padding

### Adhesive Cast Padding

**Supplier:**

Lohmann & Rauscher, Inc.

**Features & Benefits:**

- Nonwoven, polyacrylate adhesive synthetic padding
- Air permeable
- High skin tolerance
- Adhesive cast padding can be cut in strips to protect large areas, such as the trunk, pelvis or leg
- Adhesive cast edge is used around the hard edges of the cast. Preferred use for upper extremity and areas with potential pressure points.

**Specifications:**

- Materials: 50% polyester, 30% polypropylene, 20% viscose  
Adhesive = polyacrylate adhesive
- Storage: Dry area with even temperature, no higher than 25°C

**Indications:**

- Padding for immobilization after fractures, operations and orthopedic corrections

**Warranty:**

- Against manufacturers defects



Item No.	Description	Thickness	Width	Length	Unit
LR-50810-BX	Adhesive padding	5mm	19cm (7.5")	38cm (15")	10 Pcs / Bag
LR-50810-CS	Adhesive padding	5mm	19cm (7.5")	38cm (15")	60 Pcs / Case
LR-50811-BX	Adhesive padding	2mm	58cm (2')	1m (1 yd)	1 Roll / Box
LR-50811-CS	Adhesive padding	2mm	58cm (2')	1m (1 yd)	10 Rolls / Case
LR-50813-BX	Adhesive padding	2mm	28 cm (11")	10m (11 yds)	1 Roll / Bag
LR-50813-CS	Adhesive padding	2mm	28 cm (11")	10m (11 yds)	11 Rolls / Case
<b>CELLONA® CAST EDGE</b>					
LR-50812-CS	Cast edge padding	2mm	8cm (3")	2.7m (3 yds)	24 Rolls / Case