

**Ortec Adhesive**  
**Multipurpose Cement**

**Supplier:**  
Renia GmbH

**Features & Benefits:**

- Safe
- Specially designed adhesive for O&P work
- Will glue most materials used in orthotics and prosthetics
- Toluene and mek free
- Bonds Poron, EVA, plastics, polyethylenes, leather, rubber, foams, textiles, and some metals
- Dries clear
- Flexible
- Sweat and water resistant
- Can be reactivated when heated
- No separation

**Specifications:**

- Pot life: Approximately eight hours when 5–10% hardener is added
- Open assembly time: Depending on material, 10–60 minutes (TR 10–40 minutes)
- Shelf life: Minimum 12 months at 20°C

**Call Points:**

- Orthopedic shoe technician
- Pedorthist
- Foot orthotist
- Prosthetic technician
- Orthotist
- Prosthetist
- Orthotic technician



See Materials Catalog 2007/08, #26 (page 27)

1008

**Ortec Adhesive » Renia GmbH**

<b>Item No.</b>	<b>Description</b>	<b>Volume</b>
REN-170220	Contact Adhesive w/Brush	200G
REN-170264	Ortec Contact Adhesive	640G
REN-170404	Ortec Contact Adhesive	4KG (5 Litre)