


Recommended Layup for a Typical PTB

- 1 layer of Comfort Multi-Size Nylon
- 2 layers of SpectraCarb
- 1 layer of Comfort Multi-Size Nylon
- 2 layers of SpectraCarb
- 1 layer of Comfort Multi-Size Nylon
- 2 layers of SpectraCarb
- 1 layer of Comfort Multi-Size Nylon

Please keep in mind that this is a recommendation for a typical PTB layup. Adjustments in the layup will be necessary based on the weight and activity level of the patient.



Item No.	Width	Description
CSC-2	2" (5 cm)	SpectraCarb Tubular Stockinette Material
CSC-3	3" (7.5 cm)	SpectraCarb Tubular Stockinette Material
CSC-4	4" (10 cm)	SpectraCarb Tubular Stockinette Material
CSC-5	5" (12.5 cm)	SpectraCarb Tubular Stockinette Material
CSC-6	6" (15 cm)	SpectraCarb Tubular Stockinette Material
CSC-8	8" (20 cm)	SpectraCarb Tubular Stockinette Material



*Close-up view of SpectraCarb Stockinette and prosthetic.
The laminated edge can be sanded without patient nor practitioner danger of
exposure to carbon dust.*