

» TES Belt™

Transfemoral Suspension Belt

Supplier:

Syncor Ltd.

Features & Benefits:

- Features low profile, anti-roll, anti-bunch panel
- Cuff area has exposed textured neoprene on inside to securely grasp the socket
- Velcro® hook portion on one end of the waist belt is simply pressed to the exterior fabric

Specifications:

- Materials: *Inside= Neoprene with nylon*
Outside= Velcro®-compatible fabric
- Styles: *Standard and Bilateral*
- Sizes: *Eight (see chart)*
- Color: *Tan*

Indications:

- Increased suspension for transfemoral prosthetics

Call Points:

- Prosthetist
- Prosthetic technician

Warranty:

- Against manufacturers defects



TES Belt™ » Syncor Ltd.

Item No.	Description	Side	Size	Waist	Socket Brim*
SY-TES-INF-L	Transfemoral auxiliary suspension belt	Left	Infant	45–56 cm (18"–22")	23–30 cm (9"–12")
SY-TES-JUV-L	Transfemoral auxiliary suspension belt	Left	Juvenile	53–66 cm (21"–26")	28–36 cm (11"–14")
SY-TES-000-L	Transfemoral auxiliary suspension belt	Left	X-Small	68–81 cm (27"–32")	33–43 cm (13"–17")
SY-TES-100-L	Transfemoral auxiliary suspension belt	Left	Small	84–104 cm (33"–41")	41–51 cm (16"–20")
SY-TES-200-L	Transfemoral auxiliary suspension belt	Left	Medium	94–122 cm (37"–48")	48–58 cm (19"–23")
SY-TES-300-L	Transfemoral auxiliary suspension belt	Left	Large	107–135 cm (42"–53")	56–69 cm (22"–27")
SY-TES-400-L	Transfemoral auxiliary suspension belt	Left	X-Large	122–150 cm (48"–59")	66–81 cm (26"–32")
SY-TES-500-L	Transfemoral auxiliary suspension belt	Left	XX-Large	142–165 cm (56"–65")	79–91 cm (31"–36")
SY-TES-INF-R	Transfemoral auxiliary suspension belt	Right	Infant	45–56 cm (18"–22")	23–30 cm (9"–12")
SY-TES-JUV-R	Transfemoral auxiliary suspension belt	Right	Juvenile	53–66 cm (21"–26")	28–36 cm (11"–14")
SY-TES-000-R	Transfemoral auxiliary suspension belt	Right	X-Small	68–81 cm (27"–32")	33–43 cm (13"–17")
SY-TES-100-R	Transfemoral auxiliary suspension belt	Right	Small	84–104 cm (33"–41")	41–51 cm (16"–20")
SY-TES-200-R	Transfemoral auxiliary suspension belt	Right	Medium	94–122 cm (37"–48")	48–58 cm (19"–23")
SY-TES-300-R	Transfemoral auxiliary suspension belt	Right	Large	107–135 cm (42"–53")	56–69 cm (22"–27")
SY-TES-400-R	Transfemoral auxiliary suspension belt	Right	X-Large	122–150 cm (48"–59")	66–81 cm (26"–32")
SY-TES-500-R	Transfemoral auxiliary suspension belt	Right	XX-Large	142–165 cm (56"–65")	79–91 cm (31"–36")

BILATERAL

SY-TES-INF-B	Transfemoral auxiliary suspension belt	Bilateral	Infant	45–56 cm (18"–22")	23–30 cm (9"–12")
SY-TES-JUV-B	Transfemoral auxiliary suspension belt	Bilateral	Juvenile	53–66 cm (21"–26")	28–36 cm (11"–14")
SY-TES-000-B	Transfemoral auxiliary suspension belt	Bilateral	X-Small	68–81 cm (27"–32")	33–43 cm (13"–17")
SY-TES-100-B	Transfemoral auxiliary suspension belt	Bilateral	Small	84–104 cm (33"–41")	41–51 cm (16"–20")
SY-TES-200-B	Transfemoral auxiliary suspension belt	Bilateral	Medium	94–122 cm (37"–48")	48–58 cm (19"–23")
SY-TES-300-B	Transfemoral auxiliary suspension belt	Bilateral	Large	107–135 cm (42"–53")	56–69 cm (22"–27")
SY-TES-400-B	Transfemoral auxiliary suspension belt	Bilateral	X-Large	122–150 cm (48"–59")	66–81 cm (26"–32")
SY-TES-500-B	Transfemoral auxiliary suspension belt	Bilateral	XX-Large	142–165 cm (56"–65")	79–91 cm (31"–36")

* = Measure circumference of transfemoral socket around the top of the brim (as high as possible)