

**Epi-med®  
Epicondylitis Clasp**

**Vendor:**  
Thämert Ltd.

**Practitioner Benefits:**

- Easy to measure
- Easy to apply

**Patient Benefits:**

- Easy to apply
- Easy to clean

**Features:**

- Unobtrusive design does not impair bloodflow of the forearm, unlike conventional circular bandages
- Compression achieved through two pads, not a circular bandage

**Specifications:**

- The clasp can be hand washed in warm water (30°C)
- Sized by measuring the forearm approximately 4–5cm below the elbow
- Available in five sizes

**Indications:**

- Tendopathies in the elbow area caused by muscular overuse after work, sports or trauma
- Epicondylitis humeri lateralis (tennis elbow)
- Epicondylitis humeri medialis (golfer's elbow)

**Contraindications:**

- Osteochondrosis

**Call Points:**

- Occupational therapist
- Physiotherapist
- Sports medicine therapist
- Orthotist

**Warranty:**

- Against manufacture defects



**Thämert Ltd., Epi-med® Epicondylitis Clasp**

Item No.	Description	Size	Arm Circ.	Color
TH-48/001-XS	Epi-med® Epicondylitis Clasp	Extra Small	8½"–9" (22–23cm)	White w/blue trim
TH-48/001-SM	Epi-med® Epicondylitis Clasp	Small	9¼"–10" (24–25cm)	White w/blue trim
TH-48/001-MD	Epi-med® Epicondylitis Clasp	Medium	10¼"–10¾" (26–27cm)	White w/blue trim
TH-48/001-LG	Epi-med® Epicondylitis Clasp	Large	11"–11¾" (28–30cm)	White w/blue trim
TH-48/001-XL	Epi-med® Epicondylitis Clasp	Extra Large	12"–13½" (31–34cm)	White w/blue trim
TH-48/004-XS	Epi-med® Epicondylitis Clasp	Extra Small	8½"–9" (22–23cm)	Blue w/light blue trim
TH-48/004-SM	Epi-med® Epicondylitis Clasp	Small	9¼"–10" (24–25cm)	Blue w/light blue trim
TH-48/004-MD	Epi-med® Epicondylitis Clasp	Medium	10¼"–10¾" (26–27cm)	Blue w/light blue trim
TH-48/004-LG	Epi-med® Epicondylitis Clasp	Large	11"–11¾" (28–30cm)	Blue w/light blue trim
TH-48/004-XL	Epi-med® Epicondylitis Clasp	Extra Large	12"–13½" (31–34cm)	Blue w/light blue trim