

## » Alpha® Hybrid AK Liner

### Transfemoral Liner

#### Supplier:

WillowWood™

#### Features & Benefits:

- Unique gel pattern adds comfort for transfemoral amputees
- Heat moldable
- Uses Hybrid gel — a firmer durometer formulation of TPE with a Kevlar® additive
- Symmetrical gel style allows for the liner to be rotated, extending wear
- Utilizes one-way stretch fabric for controlled pistoning, as well as an accordion (ridge profile) umbrella designed to flex and conform to anatomical shapes
- Specifically designed to be used with a shuttle lock and a 10mm pin

#### Specifications:

- **Type:** *Cushion and Locking*
- **Style:** *Symmetrical*
- **Thickness:** *2.5mm (Note: Includes a 9mm distal cushion)*
- **Sizes:** *Medium Plus, Large, Large Plus and X-Large*
- **Color:** *Gray*
- New socket is required

#### Indications:


- Transfemoral applications

#### Warranty:

- 12 months, against manufacturers defects



Alpha Hybrid AK Liner  
with Symmetrical Gel Style



**Alpha® Hybrid AK Liner » WillowWood™**

TYPE	ITEM NO.	STYLE	FABRIC	COLOR	THICKNESS	SIZE
<b>CUSHION</b>	OW-H440-1494	Symmetrical	Select	Dark gray	2.5mm <sup>1</sup>	Medium Plus
	OW-H440-1496	Symmetrical	Select	Dark gray	2.5mm <sup>1</sup>	Large
	OW-H440-1498	Symmetrical	Select	Dark gray	2.5mm <sup>1</sup>	Large Plus
	OW-H440-1497	Symmetrical	Select	Dark gray	2.5mm <sup>1</sup>	X-Large
<b>LOCKING</b>	OW-H442-1494	Symmetrical	Select	Dark gray	2.5mm <sup>1</sup>	Medium Plus
	OW-H442-1496	Symmetrical	Select	Dark gray	2.5mm <sup>1</sup>	Large
	OW-H442-1498	Symmetrical	Select	Dark gray	2.5mm <sup>1</sup>	Large Plus
	OW-H442-1497	Symmetrical	Select	Dark gray	2.5mm <sup>1</sup>	X-Large

<sup>1</sup> With 9mm distal cushion

**Alpha® Hybrid AK Liner » Sizing Information**

TYPE	STYLE	SIZE	DISTAL CIRCUMFERENCE <sup>2</sup>	PROXIMAL CIRCUMFERENCE <sup>3</sup>
<b>CUSHION LOCKING</b>	Symmetrical	Medium Plus	20–32 cm (8"–12½")	30–54 cm (12"–21")
		Large	25–36 cm (10"–14")	33–64 cm (13"–25")
		Large Plus	32–47 cm (12½"–18½")	41–80 cm (16"–31½")
		X-Large	37–55 cm (14½"–22")	46–90 cm (18"–35½")

<sup>2</sup> Distal circumference is measured at 4cm (1.5") from distal end

<sup>3</sup> Proximal circumference is measured at 30cm (12") from distal end